



# High Fidelity 2X Master Mix

## Cat #: 42-501B

Contents: 500 Reactions

Storage: -20°C.

Reagent for *in vitro* laboratory use only

### General Description

**Apex** High Fidelity DNA Polymerase 2X Master Mix is an all-in-one 2x master mix containing the **Apex** High Fidelity DNA Polymerase, an optimized buffer system, dNTPs and MgCl<sub>2</sub>. Simply mix the High Fidelity DNA Polymerase 2X Master Mix with primers, DNA template and water and you are ready to carry out PCR.

**Apex** High Fidelity DNA Polymerase is a thermostable, chimeric DNA Polymerase created specifically for low-bias, high fidelity amplification of a vast range of amplicons. **Apex** High Fidelity DNA Polymerase delivers high-speed elongation and processivity, due to its fusion with a DNA-binding domain.

### Key features

- Convenient reaction setup
- High fidelity: > 60x Taq<sup>1)</sup>
- Long range amplification: 11 kb human genomic DNA
- High elongation rate: 10 sec/kb
- Excellent performance on a vast range of amplicons (High AT and high GC)
- Recommended for cloning, mutagenesis and molecular applications requiring extremely high fidelity

<sup>1)</sup> Determined through a novel NGS-based analysis of nucleotide misincorporation during PCR

### Kit Components

- Apex High Fidelity DNA Polymerase
- Optimized buffer components, 3.0 mM MgCl<sub>2</sub>
- dNTPs
- 50 mM MgCl<sub>2</sub>

### 5 M Betaine Enhancer Solution

Sold separately. Cat No.: 42-504

### Quality Control

High Fidelity DNA Polymerase is tested for contaminating activities with no traces of endonuclease activity or nicking activity. Furthermore, long range capacity is tested on a human gDNA target of 6 kb.

### Protocol

Optimal reaction conditions such as incubation times, temperatures and amount of template DNA may vary and must be determined individually. Amplification of templates with high GC content, extensive secondary structures as well as long range amplification may require more optimization - for tips see section *Strategies for Optimization*.

Prepare reaction mixtures in an area separate from that used for DNA preparation or product analysis. **Always work on ice.**

1. Thaw Apex High Fidelity 2X Master Mix and primer solutions. **It is recommended to completely thaw and thoroughly mix the master mix to ensure proper resuspension of precipitates.**
2. Prepare the reaction mix. Table 1 shows the reaction set up for a final volume of 50 µl. If desired, the reaction size may be scaled down.
3. Mix the reaction mix thoroughly and dispense appropriate volumes into reaction tubes. Mix gently, e.g., by pipetting the master mix up and down a few times.
4. Add template DNA to the individual tubes containing the reaction mix.
5. Program the thermal cycler according to Table 2. For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.
6. Place the tubes in the thermal cycler and start the reaction.

**Table 1. Reaction components for a 50 µl reaction**

Component	Vol./reaction	Final Conc.
2X Master Mix	25 µl	1X
MgCl <sub>2</sub> (50 mM)	0 – 3 µl	1.5–4.5 mM
Primer A (10 µM)	1 µl	0.2 µM
Primer B (10 µM)	1 µl	0.2 µM
Betaine (5M)*	10 - 20 µl	1 - 2M
PCR Grade Water	Variable	----
Template DNA	Variable	Variable
<b>TOTAL volume</b>	50 µl	----

\* Suggested for GC-rich amplification and long-range amplification. See section *Strategies for Optimization*.

### Strategies for PCR Optimization

#### Long-range amplification

- Longer extension times often resolve low-yield amplification of long amplicons.
- The addition of 1-2 M Betaine solution often improves reaction performance
- Increased template concentration will increase product yield.
- Increased primer concentration can increase product yield for some reactions.

### GC-rich amplification

- The addition of 1-2 M Betaine solution often improves reaction performance

### Primers

- Primers of 20 – 40 nucleotides with a GC content of 40 - 60 % are recommended.

### MgCl<sub>2</sub>

- The optimal MgCl<sub>2</sub> concentration should be determined empirically but in most cases a concentration of 1.5 mM, as provided in the High Fidelity 2X Master Mix, will produce satisfactory results. Table 3 provides the volume of 50 mM MgCl<sub>2</sub> to be added to the master mix if a higher MgCl<sub>2</sub> concentration is required.

**Table 2. Three-step PCR program for targets < 1kb**

Cycles	Duration of cycle	Temperature
1	2 min <sup>a)</sup>	98 °C
25 – 35	10 – 20 sec <sup>a)</sup>	98 °C
	15 – 30 sec <sup>b)</sup>	55 – 70 °C
	10 – 60 sec <sup>c)</sup>	72 °C
Final elongation	5 min	72 °C

<sup>a</sup> Denaturation: 2 min initial denaturation is sufficient for most templates. During thermocycling, 10 seconds usually works very well. Longer denaturation times might be required for long range PCR or amplification from templates with a high GC content.

<sup>b</sup> Primer annealing: Typically, the annealing temperature is about 3 – 5 °C below the T<sub>m</sub> (melting temperature) of the primers used. **Because of the high salt content within the High Fidelity 2X Master Mix, annealing temperature will likely be higher than with more traditional PCR buffers.**

<sup>c</sup> Extension: The recommended extension temperature is 72°C. Extension times highly depends on the length of the amplicon. **Generally, we recommend an extension time of 10-30 seconds per kb for complex genomic targets. 10 seconds per kb is often sufficient for simpler targets (such as plasmid) or short complex targets (< 3 kb). 30-60 seconds per kb is recommended for long amplicons (> 3 kb).**

**Table 3. MgCl<sub>2</sub> concentration in a 50 µL reaction**

Final MgCl <sub>2</sub> conc. in reaction (mM)	1.5	2.0	2.5	3.0	3.5	4.0	4.5
Additional volume of 50 mM MgCl <sub>2</sub>	0	0.5	1.0	1.5	2	2.5	3

### Recommended Storage and Stability

Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Optional: Store at +4 °C for up to 6 months.

### Related Products

Taq Polymerase kits (500 units)	Cat#
With 10X Standard and Ammonium Reaction Buffer	42-800B1
With 10X Combination Buffer	42-800B3
Glycerol Free	42-800B4

Hot Start DNA Polymerase kit (500 units)	Cat#
With 10X Ammonium and Combination Reaction Buffer	42-106

High Fidelity DNA Polymerase (500 units)	Cat#
With 5X High Fidelity Reaction Buffer	42-500B

Master Mixes (500 reactions)	Cat#
2X Taq RED Master Mix, 1.5 mM MgCl <sub>2</sub>	42-138
2X Taq Master Mix, Clear, 1.5 mM MgCl <sub>2</sub>	42-134
2X Hot Start Master Mix Buffer I, 1.5 mM MgCl <sub>2</sub>	42-198
2X Hot Start Master Mix Buffer I Blue, 1.5 mM MgCl <sub>2</sub>	42-144
2X High Fidelity Master Mix	42-501B

Real-time PCR (400 reactions)	Cat#
qPCR 2X Master Mix for Probe, without ROX™	42-116P
qPCR 2X Master Mix for Probe, low ROX™	42-118P
qPCR 2X Master Mix for Probe, high ROX™	42-120P
qPCR 2X GREEN Master Mix, without ROX™	42-116PG
qPCR 2X GREEN Master Mix, low ROX™	42-118PG
qPCR 2X GREEN Master Mix, high ROX™	42-120PG

Extraction Solution (500 reactions)	Cat#
DNA Extraction Solution, 500 reactions	42-503B

Genotyping PCR kit (500 reactions)	Cat#
Extract-Amp RED PCR Kit, 500 reactions	42-502B

Ultrapure dNTPs	Cat#
dNTP set, 100 mM each: 250 µl of each dA, dC, dG and dT	42-410
dNTP Mix 40 mM (1 x 500 µl): 10 mM each dA, dC, dG, dT	42-411

DNA Ladders	Cat#
Apex 100 bp-Low DNA Ladder, 250 applications	19-109
Apex 1 kb DNA Ladder, 333 applications	19-115
Apex 200 bp DNA Ladder, 200 applications	19-111
ECON Mini DNA Ladder 100-500 bp, 100 applications	19-130
ECON Low DNA Ladder 100-1000 bp, 100 applications	19-131
ECON PCR Ladder 100-3000 bp, 100 applications	19-132

Accessory reagents	Cat#
50 mM MgCl <sub>2</sub> , 3 x 1.5 ml	42-303
Nuclease-Free Water, PCR Grade, 6 x 5 ml	42-710