



Taq DNA Polymerase

Cat #: 42-801B2

Unit Size: 1000 Units

Buffers:

10X Ammonium Reaction Buffer, Mg Free
750mM Tris-HCl pH 8.5, (NH₄)₂SO₄, 1% Tween® 20.

10X Standard Reaction Buffer, Mg Free
100mM Tris-HCl, pH 8.5; 500mM KCl; 1% Tween® 20.

50 mM MgCl₂ solution

Concentration: 5 units/μl

Storage: -20°C.

Reagent for *in vitro* laboratory use only

General Description

Apex Taq DNA Polymerase is a thermostable recombinant DNA polymerase, which exhibits very high activity in primer extension and other molecular biology applications. The enzyme is isolated from *Thermus aquaticus* and has a molecular weight of approximately 94 kDa.

Apex Taq DNA Polymerase has both a 5' → 3' DNA polymerase and a 5' → 3' exonuclease activity. The enzyme lacks a 3' → 5' exonuclease activity (no proof-reading ability). Taq DNA Polymerase leaves an 'A' overhang, which makes the enzyme ideal for TA cloning.

Key Features

- The choice when high-fidelity is not required
- High performance thermostable DNA polymerase
- Ideal for rich amplifications
- Optimal for TA cloning

Unit definition

One unit is defined as the amount that incorporates 10 nmol of dNTPs into acid-precipitable form in 30 minutes at 72°C under standard assay conditions.

Storage and Dilution Buffer

Enzyme is supplied in 20mM Tris-HCl pH 8.3, 100mM KCl, 0.1mM EDTA, 1mM DTT, 0.5% Tween® 20, 50% glycerol.

Quality control

Each lot of Taq DNA Polymerase is tested for contaminating activities with no traces of endonuclease, nicking or exonuclease activity.

Protocol

This protocol serves as a guideline for primer extensions. Optimal reaction conditions such as incubation times, temperatures, and amount of template DNA may vary and must be individually determined.

- Set up reaction mixtures in an area separate from that used for DNA preparation or product analysis.

1. Thaw 10X Buffer, dNTP mix, and primer solutions. **It is important to mix the solutions completely before use to avoid localized concentrations of salts.**
2. Prepare a master mix according to Table 1. The master mix typically contains all the components needed for extension except the template DNA.

In some applications, MgCl₂ is needed for the best results. For this reason, 50mM MgCl₂ is included with the kit. Table 2 provides the volume of 50mM MgCl₂ to add to the master mix if a certain MgCl₂ concentration is required.

3. Mix the master mix thoroughly and dispense appropriate volumes into reaction tubes. Mix gently, e.g., by pipetting the master mix up and down a few times.
4. Add template DNA to the individual tubes containing the master mix.

Table 1. Reaction components (Master Mix and Template DNA) for a 50μl reaction

| Component | Vol./reaction | Final Conc. |
|----------------------------------|---------------|-----------------------|
| 10X Mg ²⁺ Free Buffer | 5μl | 1X |
| MgCl ₂ 50 mM | 0.5-5μl | 0.5-5mM |
| dNTP mix (12.5 mM of each) | 0.8μl | 0.2mM of each dNTP |
| Primer A | Variable | 0.1-1.0μM |
| Primer B | Variable | 0.1-1.0μM |
| Taq DNA Polymerase | Variable | 1 – 5 units |
| PCR Grade Water | Variable | ---- |
| Template DNA | Variable | Variable |
| TOTAL volume | 50μl | ---- |

Table 2. MgCl₂ concentration in a 50μl reaction

| Final MgCl ₂ conc. in reaction (mM) | 0.5 | 1.0 | 1.5 | 2.0 | 2.5 | 3.0 | 3.5 |
|---|-----|-----|-----|-----|-----|-----|-----|
| Additional volume of 50mM MgCl ₂ | 0.5 | 1.0 | 1.5 | 2.0 | 2.5 | 3.0 | 3.5 |


- Program the thermal cycler according to the manufacturer's instructions.

For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.

- Place the tubes in the thermal cycler and start the reaction.

Apex AGAROSE guarantee

- Low EEO (=0.12) - This means biological macro-molecules such as proteins or nucleic acids as well as larger particles such as viruses and subcellular fragments can migrate through gel's neutral properties.
- Sharp, finely resolved banding resolution as well as excellent transparency for easy reading.
- Extraordinarily low gel background after applying staining agents.
- Superior mechanical resistance for more reliable and easier handling.



| Description | Cat# |
|---------------------------------|--------|
| General Purpose Agarose (500 g) | 20-102 |

Ethidium Bromide Dropper Bottle

Adding Ethidium Bromide has never been easier and safer! For the recommended final concentration of 0.5 µg/ml, simply add one drop for every 50 ml of solution.



Concentration:
0.625mg/ml
EtBr Capacity: 5mg

| Description | Cat# |
|---------------------------|--------|
| EtBr Dropper Bottle, 10ml | 20-276 |

Ethidium Bromide Destaining Bags



Solutions and gels for safe and easy disposal. Now with our activation solution (included) our destaining bags demonstrate increased absorption efficiency speed relative to other similar products. Simply add ~5 ml of activation solution to the

| Description | Cat# |
|-------------------------------|--------|
| EtBr Destaining Bags, 25 Bags | 20-277 |

prior to use.

Related Products

| Taq Polymerase kits (500 units) | Cat# |
|--|----------|
| With 10X Standard and Ammonium Reaction Buffer | 42-800B1 |
| With 10X Combination Buffer | 42-800B3 |
| Glycerol Free | 42-800B4 |

| Hot Start DNA Polymerase (500 units) | Cat# |
|---|--------|
| With 10X Ammonium and Combination Reaction Buffer | 42-106 |

| High Fidelity - Proof reading (500 units) | Cat# |
|---|--------|
| Hi-Fi PR™ Taq 2.5 U/µl | 42-110 |

All polymerases are also available in kits, Mg²⁺ free buffers and 50 mM MgCl₂.

| Master Mixes (500 reactions) | Cat# |
|--|--------|
| 2X Taq DNA Polymerase Master Mix, 1.5 mM MgCl ₂ | 42-132 |
| 2X Taq RED Master Mix, 1.5 mM MgCl ₂ | 42-138 |
| 2X Hot Start Master Mix Buffer I, 1.5 mM MgCl ₂ | 42-198 |

The shown master mixes are ammonium based. Also available with balanced ammonium and potassium based buffers.

| Real-time PCR (400 reactions) | Cat# |
|--|----------|
| qPCR 2X Master Mix for Probe, without ROX™ | 42-116P |
| qPCR 2X Master Mix for Probe, low ROX™ | 42-118P |
| qPCR 2X Master Mix for Probe, high ROX™ | 42-120P |
| qPCR 2X GREEN Master Mix, without ROX™ | 42-116PG |
| qPCR 2X GREEN Master Mix, low ROX™ | 42-118PG |
| qPCR 2X GREEN Master Mix, high ROX™ | 42-120PG |

| Ultrapure dNTPs | Cat# |
|--|--------|
| dNTP set, 100 mM each: 250 µl of each dA, dC, dG and dT | 42-410 |
| dNTP Set, 100 mM each: 1 ml of each dA, dC, dG and dT | 42-403 |
| dNTP Mix 40 mM (1 x 500 µl): 10 mM each dA, dC, dG, dT | 42-411 |
| dNTP Mix 100 mM (2 x 1 ml): 25 mM each dA, dC, dG, dT | 42-405 |
| dNTP Mix 10 mM (10 x 1 ml): 2.5 mM each dA, dC, dG, dT | 42-406 |

Other concentrations and Single dNTPs are available.

| DNA Ladders | Cat# |
|--|--------|
| Apex 100 bp-Low DNA Ladder, 250 applications | 19-109 |
| Apex 1 kb DNA Ladder, 333 applications | 19-115 |
| Apex 200 bp DNA Ladder, 200 applications | 19-111 |
| Apex ECON Mini DNA Ladder, 100 applications | 19-130 |
| Apex ECON Low DNA Ladder, 100 applications | 19-131 |
| Apex ECON PCR Ladder, 100 applications | 19-132 |

| Accessory reagents | Cat# |
|--|--------|
| 50 mM MgCl ₂ , 3 x 1.5 ml | 42-303 |
| Nuclease-Free Water, PCR Grade, 6 x 5 ml | 42-710 |

Tween 20® is a registered trademark of ICI Americas, Inc.