

## GC10™ Competent Cells

Cat. No.	Size	Amount	pUC19 control (10 pg/μL)
42-657	0.5ml	10 x 50μl	1 x 50μl
42-658	1ml	20 x 50μl	1 x 50μl
42-659	3.75ml	75 x 50μl	1 x 50μl
42-660	1ml	5 x 200μl	1 x 50μl
42-661	15ml	75 x 200μl	1 x 50μl

Store at -70°C. FOR RESEARCH USE ONLY

### General Description

GC10 Competent Cells are chemically competent *E.coli* cells. The cells are transformed by heat-shock methods. These competent cells are resistant to the lytic bacteriophages T1 and T5.

GC10 Competent Cells are suitable for many molecular biology applications, like generating cDNA libraries from plasmid-based vectors or constructing gene banks. Blue/white screening for recombinants can be performed by including X-gal and IPTG in the agar plates.

The GC10 strain is sensitive to these common antibiotics: ampicillin, kanamycin, chloramphenicol and tetracycline. It is resistant to streptomycin.

### General Handling

- Competent cells are very sensitive to any change in temperature. Cells must be thawed on ice. The transformation should be started immediately after the cells are thawed.
- Competent cells must be treated gently. Mix cells by swirling or gently tapping the reaction tube. Do not mix by pipetting or vortexing.
- Once thawed, the cells should be used. Re-freezing thawed competent cells will result in a significant drop in transformation efficiency.

### Genotype

F- *mcrA* Δ(*mrr-hsdRMS-mcrBC*) φ80d*lacZ*ΔM15 Δ*lacX74 endA1 recA1* Δ(*ara, leu*)7697 *araD139 galU galK nupG rpsL* λ: T1R

### Efficiency

≥ 10<sup>8</sup>-10<sup>9</sup> transformants/μg of pUC19 DNA  
**pUC19 DNA Storage Buffer**

Control DNA is supplied in TE Buffer [10mM Tris-HCl pH 8.0, 0.1mM EDTA].

### Notes on Ligation Reactions

**Ligation reactions inhibit transformation.** Less transformants are observed from ligation reactions than from transformations with plasmid DNA.

Use 0.5μL of a ligation reaction per 50μL of competent cells. For best results, either purify the ligation mixture by ethanol precipitation prior to transformation or dilute the ligation reaction 3-fold in TE buffer and use 1μL per 50μL competent cells.

### Advance Preparations

- Equilibrate a non-shaking water bath to 42°C.
- Place SOC Medium at room temperature.
- Prepare LB agar plates with the appropriate antibiotic. If blue/white screening for recombinants is desired, the plates should include 40μg/mL X-gal and 1mM IPTG.
- Agar plates should be placed in a 37°C incubator for about 30 min. prior to plating.

### Transformation Protocol for Chemically Competent GC10 Cells

- Remove competent cells from -70°C and place directly in ice. Thaw cells for 5 to 10min.
- Gently mix cells by tapping tube**
- Add 1 to 50ng of DNA [or 1μL control DNA] per 50μL of cells. Swirl the pipettor tip through the cells while dispensing DNA. Gently tap tube to mix.
- Place the tubes on ice for 30min.
- Heat-shock the cells for 45sec. in a 42°C water bath. Do not shake.
- Add 450μL of room temperature SOC medium to each transformation reaction.
- Incubate at 37°C for one hour, with shaking (225 to 250rpm).
- Spread on LB agar plates containing appropriate antibiotic (e.g., 100μg/mL ampicillin for control pUC19).
- Incubate the plates at 37°C overnight (12 to 16 hours).

***SOC Medium Formulation***

2% Tryptone, 0.5% Yeast Extract, 0.4% glucose, 10mM NaCl, 2.5mM KCl, 5mM MgCl<sub>2</sub>, 5 mM MgSO<sub>4</sub>.

***Quality Control***

Cells must have a transformation efficiency of  $\geq 10^8$ - $10^9$  transformants/ $\mu$ g pUC19 DNA (non-saturating conditions).  
Cells must show resistance to T5 phage.



A Life Science Company

Corporate: ----- (800) 789.5550

Fax: ----- (888) 789.0444

e-mail: ----- support@geneseesci.com

Technical Service:----- (800) 789.5550

Web: ----- [www.geneseesci.com](http://www.geneseesci.com)