

INSTRUCTION MANUAL

Quick-DNA[™] HMW MagBead Kit Catalog No. D6060

Highlights

- High Molecular Weight DNA: Extract high molecular weight DNA up to 150 kb from any sample.
- Ultra-Pure: Highest DNA yield and purity equipped with RNA removal technology.
- Third-Generation Sequencing Ready: Optimized for long read sequencing (including Oxford Nanopore[™] and PacBio SMRT[™] sequencing).

Contents

Product Contents	1
Product Specifications	1
Product Description	2-3
Purification Guide	4
Reagent Preparation	5
Protocol	5-6
Appendix A: Enzymatic Digestion of Microbes .	7
Appendix B: Optional RNase Treatment	8
Troubleshooting Guide	9
Ordering Information	10

For Research Use Only

Ver. 1.0.0

Notes:

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product, please call 1-888-882-9682.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

[™] Trademarks of Zymo Research Corporation. Nanopore[™] is a trademark of Oxford Nanopore Technologies. Qubit[™]and NanoDrop[™] are trademarks of Thermo Fisher Scientific. TapeStation[®] is a registered trademark of Agilent Technologies, Inc.

Product Contents

Quick-DNA™ HMW MagBead Kit (Kit Size)	D6060 (96 Preps.)	Storage Temperature
Proteinase K & Storage Buffer ¹	2 x 20 mg	-20°C (after mixing)
Biofluid & Solid Tissue Buffer	25 ml	Room Temp.
Quick-DNA™ MagBinding Buffer	150 ml	Room Temp.
DNA Pre-Wash Buffer ²	50 ml	Room Temp.
g-DNA Wash Buffer	1 x 200 ml	Room Temp.
DNA Elution Buffer	50 ml	Room Temp.
MagBinding Beads	6 ml	Room Temp.

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

¹ Prior to use, reconstitute the lyophilized **Proteinase K** with 1040 μl **Proteinase K Storage Buffer**. Vortex to dissolve. Store at -20^oC.

² A precipitate may have formed in the **DNA Pre-Wash Buffer** during shipping. To completely resuspend the buffer, incubate the bottle at 30 – 37 °C for 30 minutes and mix by inversion. DO NOT MICROWAVE.

Specifications

- **Sample Types** Any cells, solid tissue, whole blood, saliva, biological fluids, buccal, swabs, stool, microbiome samples, samples stored in DNA/RNA Shield[™], *etc.*
- **DNA Purity** High quality DNA is ready for all sensitive downstream applications such as long read sequencing, PCR, endonuclease digestion, Southern blotting, genotyping, Next-Generation Sequencing, bisulfite conversion, *etc.* (*A*₂₆₀/*A*₂₃₀ ≥ 1.8).
- DNA Yield The DNA binding capacity is 10 µg per 50 µl MagBinding Beads used.
- **DNA Size** Capable of recovering genomic and mitochondrial DNA sized fragments up to 150 kb. If present, plasmid, parasitic, microbial, and viral DNA will also be recovered.
- Elution Volume 50 µl DNA Elution Buffer per 33 µl MagBinding Beads.
- Equipment: Magnetic rack, shaker and/or rotator, and microcentrifuge (optional)
- Automation: For assistance with automating/scripting this workflow onto your device, contact one of our automation experts at automation@zyomresearch.com.

Product Description

The **Quick-DNA[™] HMW MagBead Kit** is the easiest method for high molecular weight total DNA extraction (e.g., genomic, mitochondrial, viral) from any biological fluid, cell culture, solid tissue, or environmental sample. Innovative reagents and Zymo Research's unique system allow for a simple Bind, Wash, & Elute procedure that is unmatched in providing ultra-pure and high yield DNA of > 150 kb size. Isolated DNA is ready for immediate use in sensitive downstream applications including long read (Oxford Nanopore[™]) and NGS sequencing, qPCR, arrays, and methylation analysis.







High Molecular Weight DNA. Cultured *E. coli* (~10⁸ cells), cultured *L. monocytogenes* (~10⁸ cells), 50 mg human feces, and 75 µl ZymoBIOMICS[™] Microbial Community Standard (D6300) were input into the *Quick*-DNA[™] HMW MagBead kit (n=2). Length of the highest detected peak were recorded and averaged for each sample. DNA size was analyzed using Agilent's Femto Pulse system.

Notes:

For **Technical Assistance**, please contact **Zymo** at 1-888-882-9682 or E-mail tech@zymoresearch.com.



ZYMO RESEARCH CORP. Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com



Notes:

 1 If using < 200 µI sample, increase the volume to 200 µI using TE Buffer or an isotonic buffer before continuing.

¹ **DNA/RNA Shield**[™] (R1100-50) is sold separately.

Reagent Preparation

- Add 1,140 µl Proteinase K Storage Buffer to each Proteinase K (20 mg) tube prior to use. The final concentration of Proteinase K is ~20 mg/ml. Store at -20°C after mixing.
- ✓ Mix the **MagBinding Beads** until the beads are completely resuspended before use.

Protocol

The isolation consists of two steps: Sample Preparation & DNA Purification.

Sample Preparation

All steps should be performed at room temperature (20-30°C) unless specified.

Wide bore tips can be used to prevent shearing of DNA.

Biological Fluids & Cells (Whole Blood, Saliva, etc.) ≤ 200 µl

- Add 200 μl (equal volume) Biofluid & Solid Tissue Buffer to 200 μl liquid sample¹ and mix thoroughly.
- Add 20 μl Proteinase K and pipette mix 5 times. Incubate at room temperature (20-30 °C) for 20 minutes.
- 3. Proceed to DNA Purification (Page 6).

Solid Tissue (Ear/Tail Snips, Liver, Plants, *etc.*) ≤ 25 mg

- Add ≤ 25 mg solid tissue to 95 µl DNA Elution Buffer, 95 µl Biofluid & Solid Tissue Buffer, and 10 µl Proteinase K.
 Note: Plant tissue must be frozen with liquid nitrogen and ground into a fine powder before processing.
- 2. Pipette mix 5 times and incubate at 55°C for 1-3 hours or until tissue solubilizes.
- 3. Centrifuge the sample at \geq 10,000 x g with a microcentrifuge for 1 minute to pellet the debris.
- 4. Remove the supernatant while avoiding debris and transfer it to a new tube.
- 5. Proceed to DNA Purification (Page 6).

Samples in DNA/RNA Shield[™] ≤ 400 µl or 25 mg

- Add 20 μl Proteinase K to 400 μl sample in DNA/RNA Shield^{™1} and mix well. Incubate at room temperature (20-30°C) for 30 minutes. Note: Solid tissue in DNA/RNA Shield must be mechanically homogenized before processing.
- 2. Proceed to DNA Purification (Page 6).

Microbial/Environmental (Stool, Soil, Swabs) \leq 100 mg or 10⁸ Cells

1. Refer to Appendix A for enzymatic digestion of microbial samples (Page 8).

DNA Purification

- 1. Add 400 µl (equal volume) Quick-DNA[™] MagBinding Buffer to 400 µl of the sample.
- 2. Pipette mix the solution.
- 3. Add the 33 µl of MagBinding Beads¹ to each sample.
- 4. Pipette mix 5 times and place the samples on a rotator or shaker² for 10 minutes.
- 5. Transfer the sample to the magnetic stand³ until beads have separated from the solution, then remove⁴ and discard the supernatant. Transfer the sample off the magnetic stand.
- 6. Add 500 µl Quick-DNA[™] MagBinding Buffer.
- 7. Pipette mix to resuspend the beads (~5 times) and then place samples on a rotator or shaker for 5 minutes.
- 8. Transfer the sample to the magnetic stand until beads have separated from the solution, then remove⁴ and discard the supernatant. Transfer the sample off the magnetic stand.
- 9. Add 500 µl DNA Pre-Wash Buffer.
- 10. Pipette mix to resuspend the beads (~10 times).
- 11. Transfer the sample to the magnetic stand until beads have separated from the solution, then remove⁴ and discard the supernatant. Transfer the sample off the magnetic stand.
- 12. Add 900 µl g-DNA Wash Buffer and pipette mix to resuspend the beads (~10 times).
- 13. Transfer all liquid into a new microcentrifuge tube or 96-well plate.
- 14. Transfer the sample to the magnetic stand until beads have separated from the solution, and then remove⁴ and discard the supernatant. Transfer the sample off the magnetic stand.
- 15. Repeat steps 12 14.
- 16. To dry the beads, transfer the sample to a heated element and incubate for 10 minutes. If no heating element is available, air dry for 20 minutes⁵.
- 17. Add 50 µl of DNA Elution Buffer to each sample.
- Pipette mix 20 times and incubate at room temperature for 5 minutes or mix via shaker for 5 minutes at room temperature (20-30°C).
- 19. Transfer the sample to the magnetic stand until beads have separated from solution, then transfer the eluted DNA to a new tube (plate).

The eluted DNA can be used immediately or stored at \leq -20°C.

¹MagBinding Beads settle quickly, ensure that beads are kept in suspension while dispensing.

 2 Shaking speeds can be different for each shaker. Shaker should be fast enough to completely resuspending the beads (1100 – 1500 rpm).

³Magnetic stand (manual processing) or strong-field 96-well magnetic stand (*i.e.*, **ZR-96 MagStand**, P1005).

⁴ Some beads will adhere to the sides of the well. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.

⁵ Over drying the beads may result in lower DNA recovery. Beads will change in appearance from glossy black when still wet to a matte black/brown when fully dry.

 Enzymatic lysis of cells walls (e.g. Lysozyme, Zymolyase) from microbes is necessary to effectively isolate high molecular weight DNA from microbes. Fluids (Whole Blood, Saliva, Water DNA/RNA Shield, Feces) ≤ 200 μl 1. Add 100 μl (equal volume) DNA/RNA Shield[™] (2x Concentrate)¹ to up to 100 μl sample and pipette mix 10 times. Note: If sample is already resuspended in DNA/RNA Shield[™], add 100 μl DNA/RNA Shield^{™2}. 2. Incubate at room temperature (20-30°C) on a tube rotator or shaker for 5 minutes. 3. Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells 1. Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 μl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. 2. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Fluids (Whole Blood, Saliva, Water DNA/RNA Shield, Feces) ≤ 200 μl 1. Add 100 μl (equal volume) DNA/RNA Shield[™](2x Concentrate)¹ to up to 100 μl sample and pipette mix 10 times. Note: If sample is already resuspended in DNA/RNA Shield[™], add 100 μl DNA/RNA Shield^{™2}. 2. Incubate at room temperature (20-30°C) on a tube rotator or shaker for 5 minutes. 3. Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells 1. Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 μl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. 2. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Add 100 µl (equal volume) DNA/RNA Shield[™](2x Concentrate)¹ to up to 100 µl sample and pipette mix 10 times. Note: If sample is already resuspended in DNA/RNA Shield[™], add 100 µl DNA/RNA Shield^{™2}. Incubate at room temperature (20-30°C) on a tube rotator or shaker for 5 minutes. Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Note: If sample is already resuspended in DNA/RNA Shield[™], add 100 µl DNA/RNA Shield^{™2}. Incubate at room temperature (20-30°C) on a tube rotator or shaker for 5 minutes. Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Incubate at room temperature (20-30°C) on a tube rotator or shaker for 5 minutes. Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Proceed to Microbial Lysis. Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Cells and Solids (Cultured Cells, Feces, Soil, etc.) ≤ 100 mg or 10⁸ bacterial cells 1. Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. 2. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
 Resuspend up to 100 mg of sample or up to 10⁸ cells with 200 µl DNA/RNA Shield^{™2} in a microcentrifuge tube and pipette mix well. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
2. Incubate at room temperature (20-30°C) on a tube rotator for 5 minutes.
3. Proceed to Microbial Lysis.
Microbial Lysis
 Centrifuge at 5,000 x g in a microcentrifuge for 1 minute to pellet the sample. Transfer the supernatant (~180 μl) in a new microcentrifuge tube. Save both the supernatant and pellet.
Note: If a pellet does not appear after centrifugation, skip to step 9.
2. Add 100 µl PBS (user supplied) to sample pellet and pipette mix until pellet is visibly resuspended.
 Centrifuge at 5,000 x g for 1 minute to pellet the sample. Combine the supernatant with the original sample supernatant (total ~280 μl) from the previous step.
4. Add 1 ml PBS (user supplied) to the new pellet and mix until pellet is visibly resuspended.
5. Centrifuge at 5,000 x g in a microcentrifuge for 1 minute to pellet the sample and discard the supernatant.
6. Add 100 μl TE Buffer and 25 μl lysozyme ³ (100 mg/ml; user supplied) to the pellet.
7. Pipette mix until pellet is visibly resuspended, then incubate at 55°C for 30 minutes.
8. Combine the saved supernatant (~280 μ l) with the 125 μ l digested sample.
 Add 20 μl 10% SDS (user provided) and 10 μl Proteinase K. Briefly pipette mix and incubate at 55°C for 10 minutes.
10. Centrifuge 5,000 x g in a microcentrifuge for 1 minute to pellet residual debris. Transfer the supernatant to a new microcentrifuge tube.
11. Add 800 µl (2 volumes) Quick-DNA [™] MagBinding Buffer to the sample and mix well
12. Proceed to step 3 of DNA Purification on Page 6.

Appendix B: Optional RNase Treatment

This modification is necessary for downstream applications that are sensitive to trace amount of RNA. This step is not required for most long read sequencing applications (e.g., Nanopore, PacBio SMRT).

- 1. Continue after step 8 of the DNA Purification protocol on page 6.
- 2. Add 100 µl DNA Elution Buffer to the beads.
- 3. Pipette mix to resuspend the beads (~5 times).
- 4. Add 2 µl RNase A (1 mg/ml) to the sample.
- 5. Pipette mix to resuspend the beads (~5 times).
- 6. Incubate at room temperature (20 25 °C) for 20 minutes.
- 7. Add 500 µl *Quick*-DNA[™] MagBinding Buffer.
- 8. Pipette mix to resuspend the beads (~5 times) and then place samples on a rotator or shaker for 10 minutes.
- 9. Transfer the sample to the magnetic stand until beads have separated from the solution, then remove and discard the supernatant. Transfer the sample off the magnetic stand.
- 10. Continue to step 9 of the DNA Purification protocol on page 6.

Notes:

Troubleshooting Guide

Problem	Possible Causes and Suggested Solutions
Low DNA Yield or Quality	
Binding Conditions	 Binding Time. Make sure to incubate on rotator or shaker for 10 minutes after the <i>Quick-DNA</i> MagBinding Buffer has been added to the sample. Incubation for longer periods of time may help to increase yield Amount of MagBinding Beads. The volume of beads used can be increased to 50 µl and eluted in 100 µl to increase the maximum binding capacity and accommodate samples of high biomass. 33 µl is the recommended starting point and can bind up to 10 µg (sample type dependent).
Proteinase K Digestion	• Incubation Time. The optimal time is largely sample dependent. 30 minutes is recommended for liquids whereas solid tissues may be incubated overnight for complete digestion. This will maximize yields but increases protocol time.
Magnetic Beads Preparation	 Resuspension of Beads. The MagBinding Beads settle quickly. Ensure complete resuspension before use by thoroughly shaking and/or vortexing the bottle.
Poor Mixing/Resuspension	• Increase mixing cycles and/or speed. Pipette mixing of the sample is crucial for some key steps (after adding Proteinase K and after adding MagBinding Beads) to ensure sufficient resuspension. Combined tip mixing and shaking at each step is recommended for optimal DNA yields and purities.
Drying Parameters	• Prolonged Time or Increased Temperature. Over-drying beads will result in severely reduced yields. To remove residual liquid, 55°C for 10 minutes is a good starting point but can depend on specific plate dimensions and heater used.
Elution Parameters	 Low Concentration. If the final concentration of your extracted DNA is too low, use 15 µl MagBinding Beads and 30 µl DNA Elution Buffer when processing similar samples in future. Incomplete Elution. The recommended minimum elution volume is 1.5X ratio of the MagBinding Beads used (Ex. 50 µl beads to 75 µl elution). Using more volume ensures better surface coverage whereas using less volume can result in severely reduced yields and purities. Temperature Conditions. Incubating the elution step at 55°C or greater during the minutes of shaking time may increase final yield.
Low Purity	
g-DNA Wash Conditions	 New Tube Transfer. It is crucial to transfer the g-DNA Wash / MagBead mixture to a new 1.5 mL microcentrifuge tube or 96-well plate during both wash steps. This prevents salt carryover which can lower purities. Resuspension of Beads. The MagBinding Beads settle out of solution quickly, so it is important to pre-mix the beads by pipette mixing to ensure full homogeneity before additional mixing via rotator or shaker.
Elution Washes	• Insufficient Mixing. It is important to properly mix the DNA Elution Buffer when added to the MagBinding Beads. Inefficient mixing can result in lower purities.
Low Molecular Weight Size	
Rigorous Mixing	• Vortex and Shaking at High Speeds. Mixing the sample using rigorous mixing parameters (e.g., vortex and shaking at high speeds) may cause shearing of the DNA, resulting in lower size recovery. Mixing via a rotator is recommended for higher size recovery.

Ordering Information

Product Description	Kit Size	Catalog No.
Quick-DNA [™] HMW MagBead Kit	96 preps.	D6060

For Individual Sale	Amount	Catalog No.
Proteinase K & Storage Buffer	5 mg set 20 mg set	D3001-2-5 D3001-2-20
Biofluid & Solid Tissue Buffer	25 ml 100 ml	D4081-3-25 D4081-3-100
<i>Quick-</i> DNA™ MagBinding Buffer	150 ml 250 ml	D4077-1-150 D4077-1-250
DNA Pre-Wash Buffer	15 ml 30 ml 50 ml 250 ml	D3004-5-15 D3004-5-30 D3004-5-50 D3004-5-250
g-DNA Wash Buffer	50 ml 100 ml 200 ml 250 ml 400 ml	D3004-2-50 D3004-2-100 D3004-2-200 D3004-2-250 D3004-2-400
DNA Elution Buffer	1 ml 4 ml 10 ml 16 ml 50 ml	D3004-4-1 D3004-4-4 D3004-4-10 D3004-4-16 D3004-4-50
MagBinding Beads	3 ml 6 ml 8 ml 12 ml 16 ml 24 ml	D4100-2-3 D4100-2-6 D4100-2-8 D4100-2-12 D4100-2-16 D4100-2-24

Page 10

Page 11

