**General instruction for  
Micro pH electrodes  
Instrucciones generales para  
microelectrodos de pH  
Instructions g é n é rales  
pour les micro- é lectrodes de pH**

Ohaus Corporation  
7 Campus Drive  
Suite 310  
Parsippany, NJ 07054 USA  
Tel: (973) 377-9000  
Fax: (973) 944-7177

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<b>Model</b>	<b>Item NO</b>	<b>Description</b>
<b>STMICRO5</b>	30087566	Micro sample pH sensor with 5 mm diameter
<b>STMICRO8</b>	30087569	Micro sample pH sensor with 8mm diameter

**Specifications:****Connection:** BNC, Cinch**Cable length:** 1.0m**Shaft length:**

The upper part of the shaft is 50mm length with 12mm diameter; for the lower part of the shaft:

STMICRO5 is 80mm length;

STMICRO8 is 150mm length;

**Shaft diameter: (lower part of the shaft)**

5mm (STMICRO5);

8mm (STMICRO8)

**pH range:** 0-14pH**Zero potential:**  $0 \pm 20\text{mV}$ **Slope:**  $\geq 97\%$ **Impedance:**  $\leq 100\text{M}\Omega$ **Shaft material:** Glass**Reference Electrolyte:**

3M KCl solution

**Temperature range:**

0-100°C

**Reference Diaphragm:**

Sleeve or ground glass

**Micro sample pH electrodes** can be used for testing some sample in a tube, which needs the diameter of the electrode should be smaller than the tube.

**When unpacking**, check the electrode for any mechanical damage to the membrane, diaphragm or body. The length is long while the diameter is very short which make the shaft very fragile, so **be careful** when using it.


### **FAQs:**

#### **1. How do you properly store the pH electrode after every use?**

A: Rinse the electrode with pure water, put the electrode into the watering cap. If the solution in the cap is polluted or dry, clean it and fill it with 3M KCl solution.

#### **2.What's the life time?**

A: Measured samples, storage conditions and maintenance methods all could influence the life of pH electrodes; normally a pH electrode expires after 12 months.

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Ohaus Corporation, all rights reserved/ todos los derechos reservados/ tous droits réservés