

N-Gen Nitrile Exam Gloves
Chemical Resistance & Barrier Guide
Cat# 44-101

| Chemical | Nitrile |
|-------------------------------------|----------------|
| ACETALDEHYDE | P |
| ACETIC ACID | G |
| ACETONE | F |
| ACETONITRILE | F |
| ACRYLIC ACID | G |
| AMYL ALCOHOL | E |
| BENZENE | F |
| BUTYL ACRYLATE | P |
| CHROMIC ACID | F |
| DIACETONE ALCOHOL | G |
| DIESEL FUEL | E |
| DIISOBUTYL KETONE | G |
| ETHANOL | F |
| FORMALDEHYDE, 30-70% | E |
| FORMIC ACID | G |
| FREON 113 OR TF | E |
| GLUTARALDEHYDE, <5% | G |
| GLYCEROL | E |
| HEXANE | E |
| HYDROCHLORIC ACID, <30% | G |
| HYDROCHLORIC ACID, 30-70% | G |
| HYDROFLURIC ACID <50% | E |
| ISOBUTYL ALCOHOL | E |
| ISOPROPYL ALCOHOL | E |
| JET FUEL, <30% AROMATICS 73-248C | G |
| KEROSENE | E |
| METHANOL | F |
| N-PROPYL ALCOHOL | E |
| NITRIC ACID, <30% | E |
| PHENOL, >70% | G |
| PHOSPHORIC ACID, >70% | E |
| SILVER NITRATE | G |
| SODIUM CHLORIDE | E |
| SULFURIC ACID, 30-70% | F |
| XYLENES | F |

(P) Poor (F) Fair (G) Good (E) Excellent