

**E-Gen Latex Exam Gloves**  
**Chemical Resistance & Barrier Guide**  
**Cat# 43-101**

---

<b>Chemical</b>	<b>Latex</b>
ACETALDEHYDE	G
ACETIC ACID	E
ACETONE	G
ACETONITRILE	F
ACRYLIC ACID	G
AMYL ALCOHOL	G
BENZENE	F
BUTYL ACRYLATE	P
CHROMIC ACID	P
DIACETONE ALCOHOL	F
DIESEL FUEL	P
DIISOBUTYL KETONE	P
ETHANOL	E
FORMALDEHYDE, 30-70%	G
FORMIC ACID	E
FREON 113 OR TF	P
GLUTARALDEHYDE, <5%	G
GLYCEROL	E
HEXANE	P
HYDROCHLORIC ACID, <30%	E
HYDROCHLORIC ACID, 30-70%	G
HYDROFLUORIC ACID <50%	E
ISOBUTYL ALCOHOL	P
ISOPROPYL ALCOHOL	E
JET FUEL, <30% AROMATICS 73-248C	P
KEROSENE	P
METHANOL	F
N-PROPYL ALCOHOL	E
NITRIC ACID, <30%	E
PHENOL, >70%	G
PHOSPHORIC ACID, >70%	G
SILVER NITRATE	E
SODIUM CHLORIDE	E
SULFURIC ACID, 30-70%	E
XYLENES	P

---

(P) Poor      (F) Fair      (G) Good      (E) Excellent