

Formulations

PHOSPHATE BUFFERED SALINE

COMPONENTS	Cat.# 114-056 mg / L	Cat.# 114-058 Cat.# 130-058 mg / L	10X Cat.# 119-069 Cat.# 130-069 mg / L	Cat.# 114-154 Cat.# 130-154 mg / L
Inorganic Salts				
Na ₂ HPO ₄ (anhydrous)	726.00	795.00	7950.00	2160.00
KH ₂ PO ₄	210.00	144.00	1440.00	200.00
KCl				200.00
NaCl	9000.00	9000.00	90000.00	8000.00
Tween 20®				0.5 ml/L
<p>Reference: Personal Communication - Mr. George Gardiner, NIH Media Unit Tween 20® is a registered trademark of ICI Americas, Inc.</p>				