

Deep vessel with sturdy handle designed for convenient handling and pouring of liquids

- Flat bottom ensures stability and prevents tipping
- Notched pour spout design controls liquid pour
- Operating temperature range 0 to 40°C (104°F)
- Made from FDA food grade polystyrene in an ISO Class 8 clean room
- Material biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols and bases
- Recyclable, marked
- Available in 3 Sizes
 - 50mL - 100mL - 200 mL



Built-in Pour Spout
for accurate pouring

Large Sturdy Handle
for ease of movement
and pouring

200 mL



100 mL



50 mL

Flat Bottom
perfect for placing on
scales or tabletops
without tipping



Item No.	Color	Size	Volume mL	L x W x H in	L x W x H cm
120890	White	<input type="checkbox"/> Small	50	3.6 x 1.6 x 1.7	9.2 x 4.0 x 4.3
120891	White	<input type="checkbox"/> Medium	100	4.2 x 3.5 x 1.9	10.7 x 9 x 4.8
120892	White	<input type="checkbox"/> Large	200	5.8 x 4.6 x 2	14.8 x 11.7 x 5.1