


Secure & Convenient Bottle Transport

The Safety Bottle Carrier provides a safe and reliable way to transport a variety of bottled reagents in your lab. Available in 1 and 4.5 liter capacities, this carrier is the practical answer for any clinical and research lab setting

Robust Design

- Ergonomic handle prevents hand fatigue and related accidents when carrying heavier loads
- Handle swings down and away to either side for easier loading and unloading
- Centering ring stabilizes bottle in the center of carrier for safer, steadier transport
- Carriers come in four colors with designated labeling areas for easy identification
- Made of chemical resistant thermoplastic elastomer - PP blend



Movable Handle
swings down and away to allow clear access to vessel



Available in 4 Colors
for easy identification
and organization



Safety Bottle Carrier, 4.5 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120911	Red ■	7.2 x 10.5	18.3 x 26.7
120912	Black ■		
120913	Blue ■	16 H w/ handle extended	40.6 H w/ handle extended
120914	Green ■		

Safety Bottle Carrier, 1 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120915	Red ■	4.5 x 6.5	11.4 x 16.5
120916	Black ■		
120917	Blue ■	10 H w/ handle extended	25.4 H w/ handle extended
120918	Green ■		