

## Ideal for Clinical and Research laboratory staining

Tray holds 18 standard microscope slides (76 x 26 mm) securely in place for staining, rinsing and drying applications

### Robust design

- Durably constructed from solvent resistant polypropylene blend
- Non-skid Silicone feet keep tray in place

### Convenient to use

- Two ways to empty stain liquid
  - Pour spouts on each corner
  - Drain hole with plug
- Equipped with handle cut outs for easy handling

### Versatile for multiple testing formats

- Black lid included for light sensitive applications
- Humidity retention wells for specialized staining needs



Humidity retention wells



Corner pour spouts

### 18-Place Slide Staining Tray

Item No.	Color	Tray L x W x H	Cover L x W x H
120879	Black ■	16.1 x 8.9 x 1.6 in 41 x 22.5 x 4 cm	15.7 x 8.4 x 0.4 in 39.8 x 21.4 x 1.1 cm