

PRODUCT INFORMATION

Lactalbumin Hydrolysate

Cat. No. L12-121



ALPHA[™]
BIOSCIENCES

Date of Issue:
10/01/17

DESCRIPTION

A pancreatic Digest of Lactalbumin, high in concentration of amino acids, used in microbiological and tissue culture media applications.

TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS

Test Type	Specification
Dehydrated appearance	Homogeneous, free flowing powder, and yellow beige to tan
Prepared appearance	Clear to slightly hazy, gold to amber with no precipitate
Clarity	Brilliant to clear
Ash	<15%
Total Nitrogen	≥10%
pH	6.5-7.5

STORAGE

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing or if the color has changed from the original color.