

PRODUCT INFORMATION

Your Smarter Culture Media Choice



Proteose Peptone #3

Cat. No. P16-119

Date of issue
05/27/2015

DESCRIPTION

Proteose Peptone #3 is used for the cultivation of microorganisms.

Typical Physical & Chemical Characteristics:

Test Type	Specification
Appearance	tan-yellow powder
Solubility (in water)	Total at 20 g/l
Amino Nitrogen	$\geq 3.4\%$
Total Nitrogen	$\geq 10\%$
Loss on Drying	$< 6\%$
pH (2% Solution)	6.5-7.5
Ash	$< 10\%$

STORAGE

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing, or if the color has changed from the original tan-yellow.

