

Cat# 25-508B



Dulbecco's PBS

Dulbecco's phosphate-buffered saline, 1X. Without calcium and magnesium. Sterile filtered. Store at room temperature.

Components	mg/L
Potassium Chloride (KCl)	200.0000
Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄)	200.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na ₂ HPO ₄)	1150.0000