

Cat# 25-536



MEM Non-Essential Amino Acids Solution

100X solution. Sterile filtered. Store at 2° to 8°C.

Components	mg/L
Glycine (C ₂ H ₅ NO ₂)	750.0000
L-Alanine (C ₃ H ₇ NO ₂)	890.0000
L-Asparagine, Monohydrate (C ₄ H ₈ N ₂ O ₃ · H ₂ O)	1500.0000
L-Aspartic Acid (C ₄ H ₇ NO ₄)	1330.0000
L-Glutamic Acid (C ₅ H ₉ NO ₄)	1470.0000
L-Proline (C ₅ H ₉ NO ₂)	1150.0000
L-Serine (C ₃ H ₇ NO ₃)	1050.0000