



8430 Juniper Creek Lane  
San Diego, CA 92126  
[www.geneseesci.com](http://www.geneseesci.com)  
800) 789.5550

# MATERIAL SAFETY DATA SHEET

## Glycine

*In Case of Emergency: Call ChemTrac (24 hr): 800-424-9300, 703-527-3887 Call (330)425-2522*

DATE UPDATED: JANUARY 13, 2017

### SECTION 1. ----- CHEMICAL IDENTIFICATION-----

-- **Product Name** Glycine  
**Product Code(s)** 18-107, 18-108, 18-109  
**Recommended Use** For Further Manufacturing Use Only  
Not for Human or Animal Drug  
Use  
**Synonyms** AMINOACETIC ACID \* GLYCOLIXIR \* HAMPSHIRE GLYCINE

### SECTION 2. ----- HAZARDS IDENTIFICATION -----

#### Emergency Overview

##### WHMIS Classification

Not WHMIS controlled.

Not a dangerous substance according to GHS.

##### HMIS Classification

**Health hazard:** 0  
**Flammability:** 0  
**Physical hazards:** 0

##### Potential Health Effects

**Inhalation** May be harmful if inhaled. May cause respiratory tract irritation.  
**Skin** May be harmful if absorbed through skin. May cause skin irritation.  
**Eyes** May cause eye irritation.  
**Ingestion** May be harmful if swallowed.

### SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

Chemical Name	EC No.	CAS-No	Weight %
Glycine	EEC No. 200-272-2	56-40-6	95-100

### SECTION 4. ----- FIRST-AID MEASURES-----

#### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration.

#### In case of skin contact

Wash off with soap and plenty of water.

**In case of eye contact**

Flush eyes with water as a precaution.

**If swallowed**

Never give anything by mouth to an unconscious person. Rinse mouth with water.

**SECTION 5. ----- FIRE FIGHTING MEASURES -----**

**Conditions of flammability**

Not flammable or combustible.

**Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

**Special protective equipment for firefighters**

Wear self contained breathing apparatus for fire fighting if necessary.

**Hazardous combustion products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, nitrogen oxides (NOx)

**Explosion data - sensitivity to mechanical impact**

no data available

**Explosion data - sensitivity to static discharge**

no data available

**SECTION 6. ----- ACCIDENTAL RELEASE MEASURES-----**

**Personal precautions**

Avoid dust formation. Avoid breathing vapours, mist or gas.

**Environmental precautions**

Do not let product enter drains.

**Methods and materials for containment and cleaning up**

Sweep up and shovel. Keep in suitable, closed containers for disposal.

**SECTION 7. ----- HANDLING AND STORAGE-----**

**Precautions for safe handling**

Provide appropriate exhaust ventilation at places where dust is formed.

**Conditions for safe storage**

Keep container tightly closed in a dry and well-ventilated place.

**SECTION 8. ----- EXPOSURE CONTROLS/PERSONAL PROTECTION-----**

**PROTECTION Personal protective equipment**

**Respiratory protection**

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

**Hand protection**

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

**Eye protection**

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

**Skin and body protection**

Choose body protection in relation to its type, to the concentration and amount of dangerous substances,

and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

**Hygiene measures**

General industrial hygiene practice.

**Specific engineering controls**

Use mechanical exhaust or laboratory fumehood to avoid exposure.

**SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -**

**Appearance**

Form	crystalline
Colour	white

**Safety data**

pH	no data available
Melting point/freezing point	Melting point/range: 240 °C (464 °F)
Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available
Auto-ignition temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density	no data available
Water solubility	no data available
Partition coefficient: n-octanol/water	no data available
Relative vapour density	no data available
Odour	no data available
Odour Threshold	no data available
Evaporation rate	no data available

**SECTION 10. - - - - - STABILITY AND REACTIVITY - - - - -**

**Chemical stability**

Stable under recommended storage conditions.

**Possibility of hazardous reactions**

no data available

**Conditions to avoid**

no data available

**Materials to avoid**

Strong oxidizing agents

**Hazardous decomposition products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, nitrogen

oxides (NOx) Other decomposition products - no data available

**SECTION 11. ----- TOXICOLOGICAL INFORMATION -----**

**Acute toxicity**

**Oral LD50**

LD50 Oral - rat - 7,930 mg/kg

**Inhalation LC50**

no data available

**Dermal LD50**

no data available

**Other information on acute toxicity**

no data available

**Skin corrosion/irritation**

no data available

**Serious eye damage/eye irritation**

no data available

**Respiratory or skin sensitisation**

no data available

**Germ cell mutagenicity**

no data available

**Carcinogenicity**

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

**Reproductive toxicity**

no data available

**Teratogenicity**

no data available

**Specific target organ toxicity - single exposure (Globally Harmonized System)**

no data available

**Specific target organ toxicity - repeated exposure (Globally Harmonized System)**

no data available

**Aspiration hazard**

no data available

**Potential health effects**

**Inhalation**

May be harmful if inhaled. May cause respiratory tract irritation.

**Ingestion**

May be harmful if swallowed.

**Skin**

May be harmful if absorbed through skin. May cause skin irritation.

**Eyes**

May cause eye irritation.

**Signs and Symptoms of Exposure**

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**Synergistic effects**

no data available

**Additional Information**

RTECS: MB7600000

**SECTION 12. ----- ECOLOGICAL INFORMATION -----**

**Toxicity**

no data available

**Persistence and degradability**

no data available

**Bioaccumulative potential**

no data available

**Mobility in soil**

no data available

**PBT and vPvB assessment**

no data available

**Other adverse effects**

no data available

**SECTION 13. ----- DISPOSAL CONSIDERATIONS -----**

**Product**

Offer surplus and non-recyclable solutions to a licensed disposal company.

**Contaminated packaging**

Dispose of as unused product.

**SECTION 14. ----- TRANSPORT INFORMATION -----**

**DOT (US)**

Not dangerous goods

**IMDG**

Not dangerous goods

**IATA**

Not dangerous goods

**SECTION 15. ----- REGULATORY INFORMATION -----**

**WHMIS Classification**

Not WHMIS controlled.

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

**SECTION 16. ----- OTHER INFORMATION -----**

<b>Issuing Date</b>	13-Aug-2009
<b>Revision Date</b>	13-Jan-2017
<b>Revision Note</b>	No information available.
<b>Recommended Restrictions</b>	No information available

**Disclaimer**

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage,

transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

**End of MSDS**