


MATERIAL SAFETY DATA SHEET

IPTG

NFPA	HMIS	Personal Protective Equipment										
	<table border="1"> <tr><td>HEALTH</td><td>1</td></tr> <tr><td>FLAMMABILITY</td><td>1</td></tr> <tr><td>REACTIVITY</td><td>0</td></tr> <tr><td>PPE</td><td>C</td></tr> </table>	HEALTH	1	FLAMMABILITY	1	REACTIVITY	0	PPE	C			
HEALTH	1											
FLAMMABILITY	1											
REACTIVITY	0											
PPE	C											

In Case of Emergency: Call ChemTrac (24 hr): 800-424-9300, 703-527-3887 Call (330)425-2522

1. Product Identification

Synonyms: IPTG (Isopropyl-beta-D-thiogalactopyranoside) GenAR; beta-D-Galactopyranoside, 1-methylethyl 1-thio-; isopropyl beta-D-thiogalactoside
 CAS No.: 367-93-1
 Molecular Weight: 238.3
 Chemical Formula: C₉H₁₈O₅S
 Product Codes: 20-109

2. Composition/Information on Ingredients

Ingredient	CAS No	Percent	Hazardous
IPTG	367-93-1	100%	Yes

3. Hazards Identification

Emergency Overview

WARNING! TOXICOLOGICAL PROPERTIES UNKNOWN. MAY BE HARMFUL IF SWALLOWED, INHALED OR ABSORBED THROUGH SKIN. MAY CAUSE IRRITATION TO SKIN, EYES, AND RESPIRATORY TRACT

Safety Ratings (Provided here for your convenience)

Health Rating: 1 - Slight
 Flammability Rating: 1 - Slight
 Reactivity Rating: 0 - None
 Contact Rating: 2 - Moderate (Life)
 Lab Protective Equip: GOGGLES; LAB COAT; PROPER GLOVES
 Storage Color Code: Green (General Storage)

Potential Health Effects

Information on the human health effects from exposure to this substance is limited.

Inhalation:

No information found, but compound should be handled as a potential health hazard. May cause irritation to the respiratory tract. Symptoms may include coughing, sore throat, labored breathing, and chest pain.

Ingestion:

No information found, but compound should be handled as a potential health hazard. May cause irritation to the gastrointestinal tract. Symptoms may include nausea, vomiting and diarrhea.

Skin Contact:

No information found, but compound should be handled as a potential health hazard. May cause irritation with redness and pain. May be absorbed through the skin with possible systemic effects.

Eye Contact:

No information found, but compound should be handled as a potential health hazard. May cause irritation, redness and pain.

Chronic Exposure:

No information found.

Aggravation of Pre-existing Conditions:

No information found.

4. First Aid Measures

Inhalation:

Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.

Ingestion:

Give large amounts of water to drink. Never give anything by mouth to an unconscious person. Get medical attention.

Skin Contact:

Immediately flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Get medical attention. Wash clothing before reuse. Thoroughly clean shoes before reuse.

Eye Contact:

Immediately flush eyes with plenty of water for at least 15 minutes, lifting lower and upper eyelids occasionally. Get medical attention immediately.

5. Fire Fighting Measures

Fire:

As with most organic solids, fire is possible at elevated temperatures or by contact with an ignition source.

Explosion:

Fine dust dispersed in air in sufficient concentrations, and in the presence of an ignition source is a potential dust explosion hazard.

Fire Extinguishing Media:

Water spray, dry chemical, alcohol foam, or carbon dioxide.

Special Information:

In the event of a fire, wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.

6. Accidental Release Measures

Remove all sources of ignition. Ventilate area of leak or spill. Wear appropriate personal protective equipment as specified in Section 8. Spills: Clean up spills in a manner that does not disperse dust into the air. Use non-sparking tools and equipment. Reduce airborne dust and prevent scattering by moistening with water. Pick up spill for recovery or disposal and place in a closed container.

7. Handling and Storage

Freeze. Store at -15C to -21C. Isolate from incompatible substances. Containers of this material may be hazardous when empty since they retain product residues (dust, solids); observe all warnings and precautions listed for the product.

8. Exposure Controls/Personal Protection

Airborne Exposure Limits:

None established.

Ventilation System:

A system of local and/or general exhaust is recommended to keep employee exposures as low as possible. Local exhaust ventilation is generally preferred

because it can control the emissions of the contaminant at its source, preventing dispersion of it into the general work area. Please refer to the ACGIH document, Industrial Ventilation, A Manual of Recommended Practices, most recent edition, for details.

Personal Respirators (NIOSH Approved):

For conditions of use where exposure to dust or mist is apparent and engineering controls are not feasible, a particulate respirator (NIOSH type N95 or better filters) may be worn. If oil particles (e.g. lubricants, cutting fluids, glycerine, etc.) are present, use a NIOSH type R or P filter. For emergencies or instances where the exposure levels are not known, use a full-face positive-pressure, air-supplied respirator. **WARNING:** Air-purifying respirators do not protect workers in oxygen-deficient atmospheres.

Skin Protection:

Wear impervious protective clothing, including boots, gloves, lab coat, apron or coveralls, as appropriate, to prevent skin contact.

Eye Protection:

Use chemical safety goggles and/or full face shield where dusting or splashing of solutions is possible. Maintain eye wash fountain and quick-drench facilities in work area

9. Physical and Chemical Properties

Appearance: White crystalline powder.

Odor: Odorless.

Solubility: No information found.
Specific Gravity: No information found.
pH: No information found.
% Volatiles by volume @ 21C (70F): No information found.
Boiling Point: No information found.
Melting Point: 120 - 122C (248 - 252F)
Vapor Density (Air=1): No information found.
Vapor Pressure (mm Hg): No information found.
Evaporation Rate (BuAc=1): No information found.

10. Stability and Reactivity

Stability:
 Stable under ordinary conditions of use and storage.
Hazardous Decomposition Products:
 Burning may produce carbon monoxide, carbon dioxide, sulfur oxides, and nitrogen oxides.
Hazardous Polymerization:
 Will not occur.
Incompatibilities:
 Strong oxidizers.
Conditions to Avoid:
 Heat, flames, ignition sources and incompatibles.

11. Toxicological Information

For Acrylamide: oral rat LD50: 124 mg/kg; skin rabbit LD50: 1680 uL/kg; investigated as a tumorigen, mutagen, reproductive effector.

-----\Cancer Lists\-----

Ingredient	---NTP Carcinogen---		IARC Category
	Known	Anticipated	
IPTG (367-93-1) (Isopropyl-beta-D- Thiogalactopyranoside)	No	No	None

12. Ecological Information

Environmental Fate: No information found.
Environmental Toxicity:
 No information found.

13. Disposal Considerations

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations.
 Dispose of container and unused contents in accordance with federal, state and local requirements.

14. Transport Information

Not Regulated

15. Regulatory Information

-----\Chemical Inventory Status - Part 1\-----

Ingredient	TSCA	EC	Japan	Australia
IPTG (367-93-1)	Yes	Yes	No	No

-----\Chemical Inventory Status - Part 2\-----

Ingredient	--Canada--			
	Korea	DSL	NDSL	Phil.
Acrylamide (79-06-1)	Yes	Yes	No	Yes

-----\Federal, State & International Regulations - Part 1\-----

Ingredient	-SARA 302-		-SARA 313-	
	RQ	TPQ	List	Chemical Catg.
Acrylamide (79-06-1)	No	No	No	No

-----\Federal, State & International Regulations - Part 2\-----

-RCRA- -TSCA

Ingredient	CERCLA	261.33	8(d)
------------	--------	--------	------

Acrylamide (79-06-1)	No	No	No
----------------------	----	----	----

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No
SARA 311/312: Acute: Yes Chronic: No Fire: No Pressure: No

Reactivity: No (Pure / Solid)

Australian Hazchem Code: None allocated

Poison Schedule: None allocated.

WHMIS:

This MSDS has been prepared according to the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all of the information required by the CPR.

16. Other Information

NFPA Ratings: Health: 1 Flammability: 1 Reactivity: 0

Label Hazard Warning:

WARNING! TOXICOLOGICAL PROPERTIES UNKNOWN. MAY BE HARMFUL IF SWALLOWED, INHALED OR ABSORBED THROUGH SKIN. MAY CAUSE IRRITATION TO SKIN, EYES, AND RESPIRATORY TRACT.

Label Precautions:

Avoid contact with eyes, skin and clothing.

Avoid breathing dust.

Keep container closed.

Use with adequate ventilation.

Wash thoroughly after handling.

Label First Aid:

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Wash clothing before reuse. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. If swallowed, give large amounts of water to drink. Never give anything by mouth to an unconscious person. In all cases, get medical attention.

Product Use:

Laboratory Reagent.

Revision Information:

MSDS Section(s) changed since last revision of document include: 3.

Status: Accepted

Reviewed By: Brady Chung

Review Date: 7/23/2009

DISCLAIMER:

Genesee Scientific Corp. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgment in determining its appropriateness for a particular purpose. GENESEE SCIENTIFIC CORP. MAKES NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE WITH RESPECT TO THE INFORMATION SET FORTH HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS. ACCORDINGLY, GENESEE SCIENTIFIC CORP. WILL NOT BE RESPONSIBLE FOR DAMAGES RESULTING FROM USE OF OR RELIANCE UPON THIS INFORMATION.
