



# Real-Time PCR Master Mixes

## ROX Reference Dye Selection Chart

	Bio-Rad				Roche				Qiagen/Corbett		Life Technologies				Agilent Technologies			Thermo	Cepheid		
	CFX96 Touch™ & CFX384 Touch™	CFX Connect™	Opticon® 2	Chromo4™	iCycler iQ™ & MyiQ™	Lightcycler® 480	Lightcycler® 1536	Lightcycler® Nano	Lightcycler® 96	Rotor-Gene Q	Rotor-Gene 6000	7500, 7500 Fast	Vii™A7	QuantStudio™ 12K Flex	7000 7300, 7700, 7900, 7900HT	StepOne™, StepOnePlus™	Mx3000P™	Mx3005P™	Mx4000™	PikoReal™	SmartCycler®
Apex 2X GREEN Mix, High ROX™ (42-120PG, 42-121PG)															X	X					
Apex 2X GREEN Mix, Low ROX™ (42-118PG, 42-119PG)												X	X	X			X	X	X		
Apex 2X GREEN Mix, Without ROX™ (42-116PG, 42-117PG)	X	X	X	X	X	X	X	X	X	X	X			X*						X	X
Apex 2X Probe Mix, High ROX™ (42-120P, 42-121P)															X	X					
Apex 2X Probe Mix, Low ROX™ (42-118P, 42-119P)												X	X	X			X	X	X		
Apex 2X Probe Mix, Without ROX™ (42-116P, 42-117P)	X	X	X	X	X	X	X	X	X	X	X			X*						X	X

\* For openArray® experiments.