

# Gobble up these amazing BLACK FRIDAY DEALS!

**BLACK  
FRIDAY  
SALE**  
THREE WEEKS  
LONG

Disinfect tools, vessels and other research items!

## UV CLAVE™

Ultraviolet Chamber

- \* Multiple UV-C lamps and reflective walls provide 360° irradiation of chamber contents
- \* Inactivates DNA, RNA, viruses and bacteria on exposed surfaces
- \* UV-C transparent shelf eliminates blind spots
- \* Auto-shutoff when door is opened
- \* Perfect for eliminating carryover in PCR reactions



**NEW!**

31-148U	UV Clave Ultraviolet Chamber	<del>\$469.00</del> \$350.00
31-148UIC	UV-C Indicator Card, pk of 25	\$89.00
31-148UIS	UV-C Indicator Stickers, pk of 40	\$79.00

Reference Promo Code: TURKEY-RD Expires: 12/04/2020

**BLOW OUT PRICE - \$350! SAVE \$119!**



## myFUGE™ Mini

Highest Capacity Mini Centrifuge!

- \* Compare to other mini centrifuges
  - ✓ 8 x 1.5/2.0ml or 0.5ml tubes
  - ✓ 4 x PCR strips
  - ✓ 32 x 0.2ml tubes
- \* Complete with rotors/adapters for all of the above
- \* Click-n-Lok™ rotor exchange system - no tools!
- \* Improved braking, decelerates in just 1 second!

**31-500 myFuge Mini, 6000rpm/2000xg, choice of color \$289.00/2**  
Add (-B) for Blue, (-G) for Green, (-P) for Purple, (-R) for Red, (-C) for clear

**BUY ONE AT \$289.00 - GET ONE FREE!**



On-board storage  
in base for  
PCR rotor and adapters



Promo  
Code:  
TURKEY-RD

**Call Us Today!  
BLACK FRIDAY  
DEALS END  
12/4/2020**

Make the move to BLUE!


## SMARTBLUE™ Mini

Compact Blue Light Transilluminator

- \* Optimized for SmartGlow™ and SYBR™ stains
- \* Bright and uniform illumination
- \* Compact 10x10cm viewing surface
- \* Scratch resistant, glass cutting surface
- \* 2 position, 100% blue light blocking filter

31-700	SmartBlue Mini Transilluminator	\$425.00
Choose 2 below, mix or match, with 31-700 purchase - up to \$400 value		
20-278D	Apex Safe™ Loading Dye with green stain	\$200.80 FREE
20-278	Apex Safe™ PreStain with green stain	\$189.80 FREE

**2 FREE VIALS OF APEX SAFE WITH PURCHASE - UP TO A \$400 VALUE!**



Genesee Scientific Phone: 800-789-5550 email: support@geneseesci.com

# Black Friday BLOWOUT Deals



## MORTEXER™

MORE Power, MORE Flexibility... MORE Vortexer!

- \* Unique Multi-Head™ for hands free microtube vortexing
- \* Instant vortexing of tubes up to 50ml
- \* Q-Drive™ dynamic balancing system minimizes vibration
- \* Continuous or touch operation



**31-100M Mortexer Vortex Mixer w/Multi-Head \$284.00/2**

**BUY ONE AT \$284.00 - GET ONE FREE!**

Reference Promo Code: TURKEY-RD Expires: 12/04/2020

Shake things up!

## BEADBUG™

MICROTUBE HOMOGENIZER

- \* Powerful mixing action for lysis, grinding and homogenization
- \* Eliminates cross contamination
- \* Available in 3 or 6 tube capacity
- \* Cold room safe
- \* Wide variety of grinding media available



<b>31-212</b>	<b>BeadBug 3, 3x2ml tube capacity</b>	<b>\$819.00</b>
<b>31-216</b>	<b>BeadBug 6, 6x2ml tube capacity</b>	<b>\$2189.00</b>
<b>31-501</b>	<b>myFuge 12 microcentrifuge with combi-rotor</b>	<b>\$390.00</b>
<b>31-212-Z3</b>	<b>TriplePure Zirconium Beads, 3mm, pk 50 tubes</b>	<b>\$155.00</b>

Additional bead materials/sizes available



**FREE MYFUGE™ 12 WITH PURCHASE - \$390 VALUE!**

## MULTI-THERM™

COOL - HEAT - SHAKE

- \* Vortex mixer accommodates all common laboratory tubes and plates
- \* Precision contoured wells for uniform thermal transfer
- \* Exchangeable blocks for tubes 0.2ml to 50ml\*
- \* Two models: Heat/Cool or Heat only
- \* Sequence function for linking or repeating programs (HC model only)



<b>31-601</b>	<b>MultiTherm Shaker with heating/cooling</b>	<b>\$2629.00</b>
<b>31-601H</b>	<b>MultiTherm Shaker, heat only</b>	<b>\$1999.00</b>
<b>31-601CB</b>	<b>Combination block: 15x0.5ml &amp; 20x1.5ml</b>	<b>\$274.00 FREE!</b>

\*Blocks sold separately. Additional blocks available.  
See [www.benchmarkscientific.com](http://www.benchmarkscientific.com) for details



**FREE COMBINATION BLOCK WITH PURCHASE - \$274 VALUE!**

Reference Promo Code: TURKEY-RD Expires: 12/04/2020



# More Black Friday BLOWOUT Deals

**FREE  
MINIMAG  
STIRRER!**



Mix it up!

## Digital Hot Plate Stirrer

- \* Large 10x10" white ceramic surface
- \* 15L stirring capacity
- \* Digital speed and temperature control
- \* Optional probe for controlling by sample temp

31-510H 10x10" Digital Hot Plate Stirrer \$639.00  
 31-510 MiniMag Mini Stirrer \$109.00 **FREE!**  
 H3770-TP Probe for sample temp control \$109.00

Reference Promo Code: TURKEY-RD Expires: 12/04/2020



**FREE! MINIMAG STIRRER, \$109 VALUE!**

Promo  
Code:  
PM-TURKEY

## BENCHWAVER™

- \* 3D motion for superior mixing
- \* Adjustable speed (5-105rpm) and tilt (0 to ±10°)
- \* Large 13.5 x 13 inch platform
- \* Digital control

**FREE STACKING PLATFORM!**



Double your workspace!  
FREE stacking platform with purchase!

31-438 BenchWaver Nutating Mixer ~~\$1249~~ \$875.00  
 31-438-SP Stacking platform ~~\$127~~ **FREE**



**DISCOUNT PRICE AND FREE STACKING PLATFORM!**

These sales carry weight!


## Accuris Analytical and Compact Balances

### Compact Balance

- \* 500g capacity, readable to 0.01g
- \* Perfect for general weighing applications

### Analytical Balance

- \* 210g capacity model, readable to 0.0001g
- \* 120g capacity model, readable to 0.00001g
- \* Touch screen (120g) or LCD (210g)
- \* RS232 for data output



31-455 Accuris Compact Balance, 500g ~~\$305.00~~ \$240.00  
 31-444A Accuris Analytical Balance, 210g ~~\$2030.00~~ \$1050.00  
 W3002A-120 Analytical Balance Tx, internal calibration, 120g \$3175.00 **PLUS 2 FREE Compact balances**

**Compact - Regular price \$305, BLOWOUT Price \$240!**  
**Analytical - Regular price \$2030.00, BLOWOUT Price \$1050!**  
**Analytical Tx - 2 FREE Compact balances with purchase!**

Free offers require purchase of designated items at listed price. Offer expires 12/04/20. Not valid on previous purchases. Cannot be combined with other promotions. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY-RD.

**Genesee Scientific Phone: 800-789-5550 email: support@geneseesci.com**