

SARTORIUS

Vivaspin® Turbo 15 PES

Purchase 2 Vivaspin® Turbo 15 PES , Get a 40% Discount

Enabling the Fastest Sample Concentration with Highest Recoveries

Vivaspin® Turbo 15 PES and RC centrifugal concentrators offer the optimal solution to any concentration or buffer exchange application with their broad range of MWCOs and unique membrane choice. Highest flow rates are achieved due to their twin vertical membranes which minimize protein polarization and subsequent fouling of the membrane. Additionally, their sleek internal profile ensures maximum process speeds right down to the last 100 μ L.

Benefits:

- Complete Recoveries
- Fastest Concentrations
- Highest Chemical Resistance
- Unique PES and RC Device Range
- Optimal Sample to Membrane Compatibility



Get a 40% Discount

Vivaspin® Turbo 15 PES

Purchase 2 Vivaspin® Turbo 15 PES , Get a 40% Discount



Purchase 2 Vivaspin® Turbo 15 PES , Get a 40% Discount

Purchase any two Vivaspin Turbo products listed below and receive a 40% discount

Model Number	Genesee Part no.	Description	List Price	Special Pricing
VS15T01	60-965	Vivaspin® Turbo 15 Centrifugal Concentrator Polyethersulfone, 12 pc treatment Cartridge, 10kDa	\$191.00	\$114.60
VS15T11	60-967	Vivaspin® Turbo 15 Centrifugal Concentrator Polyethersulfone, 12 pc treatment Cartridge, 5 kDa	\$191.00	\$114.60
VS15T21	60-968	Vivaspin® Turbo 15 Centrifugal Concentrator Polyethersulfone, 12 pc treatment Cartridge, 30 kDa	\$191.00	\$114.60
VS15T31	60-970	Vivaspin® Turbo 15 Centrifugal Concentrator Polyethersulfone, 12 pc treatment Cartridge, 50 kDa	\$191.00	\$114.60
VS15T41	60-972	Vivaspin® Turbo 15 Centrifugal Concentrator Polyethersulfone, 12 pc treatment Cartridge, 100 kDa	\$191.00	\$114.60

Special offer valid only for the US until December 31, 2024. Cannot be combined with any other offer.

Contact:

