

SARTORIUS

Special Entris® II Offer

Best Value in its Class for Basic Weighing Tasks

Every Sartorius balance offers quality, value, and consistency

No matter what you're weighing, the new Entris® II balance is always the right choice. Offering unrivalled value and backed by more than 150 years of German engineering expertise, the Entris® II comes in two product lines, so you can find the balance that meets your specific weighing needs.

Our Entris® II Advanced line gives you additional value to the Entris® II Essential portfolio, with over 38 models offering benefits such as real-time level support, integrated protection systems, CalAuditTrail, a graphic touch display and 13 built-in application programs.



Save up to 30%

**Genesee
Scientific**
A Life Science Company

Entris® II Laboratory Balance

Main Benefits & Features

- Fastest stabilization time in its class – using state- of-the-art weighing sensors
- 12 built-in applications (13 for the Advanced line) with GLP|GMP compliant printout|data output
- Polybutylene Terephthalate (PBT) housing for high chemical resistance
- Plug-and-play functionality for automatic connection to a PC
- Built-in real-time level support and integrated protection systems (only Advanced line)



Special Entris® II Offer Up to 30% off select Balances

Entris® II Advanced Balance

Model Number	Genesee Part no.	Description	List Price	Special Pricing
BCA2201-1S	60-177	Entris II Advanced Balance 2200g 0.1g	\$1,660.00	\$1,162.00
BCA2201i-1S	60-178	Entris II Advanced Balance int.cal. 2200g 0.1g	\$1,990.00	\$1,393.00
BCA223-1S	60-181	Entris II Advanced Balance 220g 1mg	\$1,990.00	\$1,393.00
BCA223i-1S	60-182	Entris II Advanced Balance int.cal. 220g 1mg	\$2,330.00	\$1,631.00
BCA224-1S	60-183	Entris II Advanced Balance 220g 0.1mg	\$4,090.00	\$2,863.00
BCA224i-1S	60-184	Entris II Advanced Balance int.cal. 220g 0.1mg	\$4,420.00	\$3,094.00
BCA3202-1S	60-185	Entris II Advanced Balance 3200g 10mg	\$2,430.00	\$1,701.00
BCA3202i-1S	60-186	Entris II Advanced Balance int.cal. 3200g 10mg	\$2,770.00	\$1,939.00
BCA323-1S	60-187	Entris II Advanced Balance 320g 1mg	\$2,220.00	\$1,554.00
BCA323i-1S	60-188	Entris II Advanced Balance int.cal. 320g 1mg	\$2,540.00	\$1,778.00
BCA324i-1S	60-189	Entris II Advanced Balance int.cal. 320g 0.1mg	\$4,860.00	\$3,402.00
BCA423-1S	60-192	Entris II Advanced Balance 420g 1mg	\$2,430.00	\$1,701.00
BCA423i-1S	60-193	Entris II Advanced Balance int.cal. 420g 1mg	\$2,770.00	\$1,939.00
BCA5201-1S	60-194	Entris II Advanced Balance 5200g 0.1g	\$1,880.00	\$1,316.00
BCA5201i-1S	60-195	Entris II Advanced Balance int.cal. 5200g 0.1g	\$2,220.00	\$1,554.00
BCA8201-1S	60-202	Entris II Advanced Balance 8200g 0.1g	\$2,330.00	\$1,631.00
BCA8201i-1S	60-203	Entris II Advanced Balance int.cal. 8200g 0.1g	\$2,650.00	\$1,855.00

Entris® II Essential Balance

Model Number	Genesee Part no.	Description	List Price	Special Pricing
BCE124-1S	60-210	Entris II Essential Balance 120g, 0.1mg	\$2,770.00	\$1,939.00
BCE124i-1S	60-211	Entris II Essential Balance int cal 120g, 0.1mg	\$2,990.00	\$2,093.00
BCE2201-1S	60-212	Entris II Essential Balance 2200g, 0.1g	\$1,390.00	\$973.00
BCE2201i-1S	60-213	Entris II Essential Balance intcal 2200g, 0.1g	\$1,670.00	\$1,169.00
BCE2202-1S	60-214	Entris II Essential Balance 2200g, 10mg	\$1,630.00	\$1,141.00
BCE2202i-1S	60-215	Entris II Essential Balance int cal 2200g, 10mg	\$1,920.00	\$1,344.00
BCE224-1S	60-219	Entris II Essential Balance 220g, 0.1mg	\$3,130.00	\$2,191.00
BCE224i-1S	60-220	Entris II Essential Balance int cal 220 0.1mg	\$3,360.00	\$2,352.00
BCE323-1S	60-225	Entris II Essential Balance 320g, 1mg	\$1,630.00	\$1,141.00
BCE323i-1S	60-226	Entris II Essential Balance int cal 320g, 1mg	\$1,920.00	\$1,344.00
BCE4202-1S	60-228	Entris II Essential Balance 4200g, 10mg	\$2,010.00	\$1,407.00
BCE4202i-1S	60-229	Entris II Essential Balance int cal 4200g, 10mg	\$2,280.00	\$1,596.00
BCE6202-1S	60-238	Entris II Essential Balance 6200g, 10mg	\$2,640.00	\$1,848.00
BCE6202i-1S	60-239	Entris II Essential Balance int cal 6200g, 10mg	\$2,880.00	\$2,016.00
BCE64-1S	60-246	Entris II Essential Balance 60g, 0.1mg	\$2,390.00	\$1,673.00
BCE64i-1S	60-247	Entris II Balance int cal 60g, 0.1mg	\$2,640.00	\$1,848.00

Special offer valid only for the US to December 31st, 2024.

Contact: