

Trade-in Your Old Rotors for a New Beckman Coulter Life Sciences Rotor

Offer Ends
12.31.2024

Contact your sales representative to calculate your savings on any Beckman Coulter Life Sciences Rotor or purchase a J-Lite rotor and get a HarvestLine System at no charge (see back).

Even the most well-designed and maintained rotor will eventually reach the point where it could fail during operation. Corrosion, stress corrosion and metal fatigue—invisible to the naked eye—eventually take their toll. Consider retiring any rotor that has reached the end of its warranty period. Even if damage is contained, a rotor failure can mean the loss of valuable samples and severely damaged instruments—a gamble that could be costly.

Schedule your rotor safety training today!

High g-Force

JS-24.38
Swinging-Bucket
Max RCF:
103,860 x g**



Capacity

J-Lite JLA-8.1000
Fixed-Angle
Max Capacity: 6 L[†]



Performance

JA-14.50 Fixed-Angle
Max Capacity:
16 x 50 mL conical tubes[†]
(Dual-locking version is
BioCertified.)



Versatility

JS-5.3 Swinging-Bucket
Max Capacity: 24
microplates, 8 deep-well
plates, 500 mL conical
bottles, (132) 13 mm dia.
tubes[†]



(BioCertified when used in
conjunction with bucket covers
or aerosol covers.)

- 30 rotor options**
- Compatible with existing Beckman Coulter High Performance Centrifuges**
- Swinging-Bucket and Fixed-angle rotors capable of greater than 100,000 x g***

- BioSafe & Dual Locking options
- Reduce Spin Time with improved k factor due to Efficient Engineering ¥¥
- J-Lite weight rotor options**

ACCELERATING
answers



Light on Weight, Heavy on Performance

Stop compromising the integrity of your sample for a “lighter” rotor. Beckman Coulter J-Lite rotors combine lightweight design, large capacity, performance and safety.

Benefits of J-Lite Rotors

- Lightweight and easy to install, reducing strain on laboratory personnel
- Improved k factor with shorter and more efficient runs
- Enhanced biosafety to minimize lab personnel exposure to biohazards
- Optimize harvesting workflows with Beckman Coulter HarvestLine Bottle Liner

PROMOTION – Purchase a new J-Lite rotor and get a HarvestLine System Liner Kit at no charge***

Watch the HarvestLine System Video

***Offer only valid for rotors that use 1 liter bottles



Contact your Life Sciences Sales Representative for details:

Name:

Email:

Phone:

* When used in the Avanti JXN-30.

** Please refer to the centrifuge user manual for the list of compatible rotors.

*** JS-24.38, JS-24.15 and JA-30.50 Ti rotors with the Avanti JXN-30.

BioCertified is a term used to describe our products which have been tested and validated to demonstrate containment of microbiological aerosols by an independent, third-party facility (Health Protection Agency, Porton Down, UK or USAMRIID, Ft. Detrick, MD, USA). Improper use or maintenance may affect seal integrity and thus, containment.

BioSafe and BioSafety are terms intended to describe the enhanced biocontainment features of our products.

¥ Please refer to the rotor user manual for the list of tubes and adapters.

¥¥ Please refer to the rotor IFU for K-factor information.



© 2024 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries. All other trademarks are the property of their respective owners.

For Beckman Coulter's worldwide office locations and phone numbers, please visit Contact Us at beckman.com
2024-AMER-EN-105146-V2 (V2.0)

