

GenClone® Cell Strainers

Pore Sizes Ranging from 5µm to 100µm

- 8 pore sizes for greater cell sorting precision
- Advanced design - faster flow through
- Included adapter for tubes smaller than 50mL
- Compatible with 6-well plates
- Interlocking and stackable
- 1.5mm feet for use in dishes & plates



5µm Strainer
25-381 - Orange



10µm Strainer
25-382 - Teal



20µm Strainer
25-383 - Red



30µm Strainer
25-384 - Purple



40µm Strainer
25-375 - Blue



50µm Strainer
25-385 - Grey



70µm Strainer
25-376 - White



100µm Strainer
25-377 - Yellow

Cat No.	Description	Qty
GenClone® Advanced Cell Strainers		
25-381	Cell Strainer, 5µm, Orange, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-382	Cell Strainer, 10µm, Teal, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-383	Cell Strainer, 20µm, Red, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-384	Cell Strainer, 30µm, Purple, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-375	Cell Strainer, 40µm, Blue, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-385	Cell Strainer, 50µm, Grey, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-376	Cell Strainer, 70µm, White, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-377	Cell Strainer, 100µm, Yellow, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
GenClone® Cell Strainer Accessories		
25-278	Reducing Adapter, Sterile, Fits 5mL, 15mL, 12×75mm, and 17×100mm	25/pk
25-280	PolyPestle™ Double-Sided Pestle for Cell Strainers	50/pk
21-395	Centrifuge Tube, Attached Cap, Sterile, Self-Standing, Bags of 25	500/cs

Reducing Adapter - 25-278



Allows coupling with 5mL, 15mL,
12×75mm and 17×100mm Tubes

Individual Blister Packaging



Sterile; RNase, DNase, pyrogen
and cytotoxin free



Poseidon™ Digital Electronic Pipette Controller

- Enhanced flow control – Max aspirate and dispense rates
- Independent aspiration and dispense speed settings
- Pressure sensitive triggers for even more control
- Improved stand - Locks in place with a click
- USB-C charge port – No searching for specific chargers



Digital Controller
33-901D

Pipette Controllers

Manual, Electronic, Long-Neck

Manual Controller 33-902

No batteries required
Autoclavable lower body
Single-lever controls

Electronic Controller 33-901

Retractable kickstand
Contoured handle
Available w/ US or EU plug

Long-Neck Controller 33-901XL

Reduced arm strain
Perfect for under the hood
Available w/ US or EU plug



Cat No.	Description	Qty
Poseidon Digital Pipette Controller		
33-901D	Poseidon™ Digital Pipette Controller, Retractable Stand, USB-C Plug	1 each
Other Poseidon Pipette Controllers		
33-901	Poseidon™ Pipette Controller, Retractable Stand, US Type Plug	1 each
33-901XL	Poseidon™ Pipette Controller XL, Retractable Stand, US Type Plug	1 each
33-902	Poseidon™ Manual Pipette Controller, Single Lever Controller	1 each



BenchBin™ Biohazard Bin

Benchtop Wastebin with Spring Lid

- Safely dispose of tips, tubes, gloves, etc.
- Secures and covers hazardous waste
- Easy push open/close spring lid
- Compliant with OSHA standard 29 CFR 1910.1030
- Compact footprint – Minimize space on the bench



Cat. No.	Description	Qty
30-165	BenchBin™, 7×5×7", 1.5L Capacity, fits 30-162R Autoclave Bags	1 each
30-166	BenchBin 3L™, 6×5×8½", 3L Capacity, fits 30-166R Autoclave Bags	1 each
30-162R	Autoclave Bags, 8½×11" Printed Markings, fits 30-165 BenchBin™	100 / pack
30-166R	Autoclave Bags, 10×12½" Printed Markings, fits 30-166 BenchBin 3L™	100 / pack
30-165K	Starter Kit: 1 BenchBin™ 1.5L (30-165) and 400 Autoclave Bags (30-162R)	1 each
30-166K	Starter Kit: 1 BenchBin 3L™ (30-166) and 400 Autoclave Bags (30-166R)	1 each

BowTie™ Biohazard Bags & Bins

Drawstring Bags and Hands-Free Bins



- Easy drawstring closure
No wire tie needed!
- 12 or 30 gallon capacity bags
- Color-changing sterilization indicator
- Hands-free foot pedal lid on bins



Cat. No.	Description	Qty
BowTie™ Biohazard Bags - Drawstring, Marking Area, Sterilization Indicator		
30-210	27×35", Fits 30-211 and 30-219 Biohazard Bin	100 / pack
30-218	33×41", Fits 08-026 Biohazard Bin	100 / pack
BowTie™ Biohazard Bins - Attached Lid and Hands-Free Foot Pedal		
30-211	8 Gallon, 10×14×19", Fits 30-210 Biohazard Bag	1 Each
30-219	13 Gallon, 16×16×24", Fits 30-210 Biohazard Bag	1 Each
08-026	26 Gallon, 19×19×33", Fits 30-218 Biohazard Bag	1 Each

Olympus Plastics™ FACS Tubes

Tubes and Strainer Caps for Flow Cytometry

Material

Tube: Polystyrene
Cap Frame: HDPE
Cap Mesh: Nylon

- Compatible with FACS™ (BD) instruments
- Strainer Caps: 20µm, 35µm, 70µm, 100µm mesh
Standard, non-strainer caps also available
- 12×75mm (5mL) tubes in crystal clear polystyrene



Cat No.	Description	Qty
Strainer and Standard Caps without Tubes		
28-157	20µm, Red, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-155	35µm, Green, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-158	70µm, White, Strainer Cap, Opaque, Sterile, Bags of 25	500 / pack
28-159	100µm, Yellow, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-152	Standard Cap, Non-Strainer, Non-Sterile, Bags of 125	1000 / case
FlowTubes™ - FACS Tubes with Included Caps		
28-154	5mL, 12×75mm, 35µm Strainer Cap, Sterile, Bags of 25	500 / case
28-156	5mL, 12×75mm, 35µm Strainer Cap, Sterile, Individual Wrap	500 / case
28-150	5mL, 12×75mm, Standard Cap, Sterile, Bags of 25	500 / case
FlowTubes™ - FACS Tubes without Caps		
28-151	5mL, 12×75mm, Without Cap, Sterile, Tray Packs of 250	1000 / case
28-153	5mL, 12×75mm, Without Cap, Non-Sterile, Bulk Econo-Pack	4000 / case




Centrifuge Tubes

5mL to 500mL Capacities

- Perfect size for any application
- RNase, DNase, & pyrogen free
- Vivid graduations

Temperature
Min: -86°C
Max: 121°C

Material
Tube: Polypropylene
Cap: Polyethylene



Cat No.	Description	Qty
5mL MacroTubes®		
24-285SC	Separate Caps, Non-Sterile, Bulk Bagged	500 / case
24-285SCB	Attached Caps, Sterile, Bags of 100	500 / case
24-285SCS	Attached Caps, Sterile, Racks of 50	500 / case
24-285SBK	Attached Caps, Opaque Black, Sterile, Racks of 50	500 / case
24-285V	Attached Caps, Vivid Graduations, Sterile, Racks of 50	500 / case

10mL Centrifuge Tubes		
21-412	Attached Cap, Sterile, Bags of 50	1000 / case

15mL Centrifuge Tubes		
21-408	Attached Cap, Sterile, Racks of 50	500 / case
21-408B	Attached Cap, Sterile, Bags of 50	500 / case
21-419	Attached Cap, Sterile, Eco-Plastic Racks of 25	500 / case
21-400	Attached Cap, Sterile, Individually wrapped	300 / case

Cat No.	Description	Qty
25mL Centrifuge Tubes		
24-295	Attached Cap, Sterile, Racks of 25	200 / pack
24-296	Attached Cap, Sterile, Bags of 25	200 / pack

30mL Centrifuge Tubes		
21-402	Screw Cap, Self-Standing, Sterile, Bags of 25	500 / case

50mL Centrifuge Tubes		
21-409	Attached Cap, Sterile, Racks of 25	500 / case
21-409B	Attached Cap, Sterile, Bags of 25	500 / case
21-421	Separate Screw Cap, Non-Sterile, Bulk Bag	500 / case
21-420	Attached Cap, Sterile, Eco-Plastic Racks of 25	500 / case
21-401	Attached Cap, Sterile, Individually Wrapped	300 / case
21-395	Attached Cap, Sterile, Self-Standing, Bags of 25	500 / case

Cat No.	Description	Qty
Large Volume Centrifuge Tubes		
21-422	175mL, Sterile PP, Conical Bottom, Bags of 6	120 / case
21-423	225mL, Sterile, PP, Conical Bottom, Bags of 5	90 / case
21-404	250mL, Sterile, PP, Conical Bottom, Bags of 5	70 / case
08-068	500mL, Sterile, PP, Conical Bottom, Bags of 5	40 / case

Cat No.	Description	Qty
Centrifuge Tube Racks		
27-470	4x50/25mL, 12x15/10/5mL, 32x1.5/0.5mL Tubes	5 / pack
27-117	Decanting, 50x5mL, Silicone Grip Wells	1 each
27-386	30 Slots for 10/15mL, 20 Slots for 25/50mL	5 / pack
27-471	6 Slots for 175, 225, or 250mL Tubes	5 / pack
24-323	4 Slots for 500mL Tubes	1 each

TotalPure™ Tubes

Metal-Free Centrifuge Tubes

- Ultra-Pure PP certified free of trace metals
- Hi-Viz™ graduations with max freeze line
- Sterile - resealable bags and individually wrapped



Cat No.	Description	Qty
21-425	15mL, Metal Free, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
21-426	50mL, Metal Free, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
21-427	15mL, Metal Free, Sterile, Hi-Viz™ Graduations, Individually Wrapped	500 / case
21-428	50mL, Metal Free, Sterile, Hi-Viz™ Graduations, Individually Wrapped	500 / case

ProSeries™ Tubes

Performance Centrifuge Tubes

- 20,000xg Max RCF marked directly on tube
- FreezeTech™ indicator changes color at 0°C
- Hi-Viz™ graduations with max freeze line
- Resealable sterile bags & eco-plastic racks

Eco-Plastic Racks



Autoclavable, recyclable,
easily flattened after use

**BUY 2
GET 1
FREE**

20,000xg Max RCF
**20K
MAX x g**

FreezeTech™
Temperature Indicator
**BLUE VISIBLE
AT 0°C**

Cat No.	Description	Qty
21-415	15mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500 / case
21-416	15mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
21-417	50mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500 / case
21-418	50mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case

Promo Code: RE2025

Disposable Cuvettes

Spectrophotometry Cuvettes

- Fits most spectrophotometers and photometers
- Polystyrene for use in visible range (340-800nm)
- Protective packaging – dust-proof and non-scratch
- Absorption reproducibility error of less than $\pm 1\%$



Cat No.	Description	Qty
21-190	Cuvette, PS, 1.5mL semi-micro, round opening, 10mm light path	1000 / pack
21-191	Cuvette, PS, 2.5mL semi-micro, round opening, 10mm light path	1000 / pack
21-192	Cuvette, PS, 1.5mL semi-micro, square opening, 10mm light path	1000 / pack
21-193	Cuvette, PS, 4.5mL macro, square opening, 10mm light path	1000 / pack
21-194	Cuvette Plug Cap, Round Closure, fits S4100-1 and S4100-2 cuvettes	500 / pack
21-195	Cuvette Plug Cap, Square Closure, fits S4100-3 and S4100-4	500 / pack



Solution Reservoirs

Standard 12-Channel Width



- Graduated inner walls
- Micro-Trough U-shaped bottoms
- Pouring spouts on all 4 corners
- RNase, DNase, & pyrogen free
- Pack of 5 protective lids with standard 25 and 50mL reservoirs


Description	25mL Standard*	50mL Standard	50mL Divided	100mL Standard	QTY
Non-Sterile Bulk Pack	28-134	28-109	28-160	28-163	100 / case
Sterile Individually Wrapped	28-134B	28-110	28-161	28-164	100 / case
Sterile Packs of 5	28-134C	28-111	28-162	28-165	200 / case
Lids - Sterile Packs of 5	28-135	28-111LID	N/A	N/A	200 / case

*5 FREE Lids included with every purchase of Standard 25mL and Standard 50mL reservoirs.

ClearSeal™

Screw Cap Microtubes

- Leak-proof O-Ring seal
- Ultra-clear polypropylene
- Resealable sterile bags of 50
- RNase, DNase, and Pyrogen free
- White printed graduations
- Conical or standing bottoms



Temperature	Material	Max RCF
Min: -86°C Max: 121°C	Polypropylene	22,000xg

Cat No.	Description	Qty
24-292	0.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000 / case
24-292B	0.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap,	1000 / case
24-297	1.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000 / case
24-297B	1.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000 / case
24-298	2.0mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap, Graduated	1000 / case
24-298B	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000 / case
24-299	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Loop Screw Cap, Graduated	1000 / case

To receive promotional packages, please reference **Promo Code RE2025** on your order. All promotions are valid while supplies last and expire **December 31, 2025**. All prices, specifications, and offers effective **April 1, 2025**. Subject to change without notice.

OptiWell™ Consumables

For KingFisher™ Platforms



Deep Well Plates

Printed Alphanumerics
Virgin, medical grade PP
24 or 96 square well models

Sealing Mats

Leak-proof silicone sealing
Printed Alphanumerics
Minimize sample loss

Tip Combs

Medical-grade PP
KingFisher™ Flex compatible
8, 24 or 96 tip models

Cat No.	Description	Qty
Deep Well Plate - Square Well, V-Bottom, PP, Printed Alpha-numerics		
30-401	Deep Well Plate, 24 Square Well, 10mL, 41.9mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
30-415	Elution Plate, 96 Square Well, 0.5mL, 15.1mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
30-425	Deep Well Plate, 96 Square Well, 1.0mL, 24.4mm Height, V-Bottom, Sterile, for KingFisher™	100 / case
30-445	Deep Well Plate, 96 Square Well, 2.2mL, 44.2mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
Sealing Mat - Square Well, Silicone, Printed Alpha-numeric		
30-467	Sealing Mat, 24 Square Well, Silicone, Sterile	50 / case
30-455	Sealing Mat, 96 Square Well, Ø8.4mm, "+" Opening, Alphanumerics, Silicone, Sterile	50 / case
30-457	Sealing Mat, 96 Square Well, Ø8.4mm, "-" Opening, Alphanumerics, Silicone, Sterile	50 / case
Tip Comb - V-Bottom, PP		
30-473	Tip Comb, 8 Tip, V-Bottom, Sterile, PP, for KingFisher™	480 / case
30-475	Tip Comb, 24 Tip, V-Bottom, Sterile, PP, for KingFisher™	50 / case
30-477	Tip Comb, 96 Tip, V-Bottom, Sterile, PP, for KingFisher™	100 / case

Visit geneseesci.com for non-sterile models and dozens of other plate/mat configurations.

KingFisher™ is a registered trademark of Thermo Fisher Scientific

Olympus Reservoirs

Reduced 8-Channel Width

- Specifically designed for 8-channel pipettes
- 75mm inner chamber
- Conserve valuable reagents



Description	5mL 8-Channel	10mL 8-Channel	25mL 8-Channel	QTY
Non-Sterile Bulk Pack	28-112	28-115	28-131*	300 / case*
Sterile Individually Wrapped	28-113	28-116	28-132	100 / case
Sterile Packs of 5	28-114	28-117	28-133	200 / case

*28-131 has quantity of 100/pack

PCR Sealing Films & Foils

Sterile Seals for PCR, qPCR, ELISA, EIA

- All types available sterile or non-sterile
- Evaporation rates < 5%
- 10 inner bags of 10, 100 per pack
- RNase, DNase, DNA, & PCR inhibitor free



Cat No.	Description	Qty
12-205	PCR Sealing Foil, Aluminum Foil, Adhesive, Pierceable	100 / pack
12-205S	PCR Sealing Foil, Aluminum Foil, Adhesive, Pierceable, Sterile	100 / pack
12-207	Sealing Film, Transparent, Easy Peel Low Adhesive, for ELISA & EIA	100 / pack
12-207S	Sealing Film, Transparent, Easy Peel Low Adhesive, for ELISA & EIA, Sterile	100 / pack
12-209	PCR Sealing Film, Transparent, Adhesive	100 / pack
12-209S	PCR Sealing Film, Transparent, Adhesive, Sterile	100 / pack
12-210	qPCR Sealing Film, Optically Clear, Standard Adhesive	100 / pack
12-210S	qPCR Sealing Film, Optically Clear, Standard Adhesive, Sterile	100 / pack
12-212	qPCR Sealing Film, Optically Clear, Non-Tacky Pressure Activated Adhesive	100 / pack
12-212S	qPCR Sealing Film, Optically Clear, Non-Tacky Pressure Activated Adhesive, Sterile	100 / pack