


## ProSeries™ Tubes

### Performance Centrifuge Tubes

- 20,000xg Max RCF marked directly on tube
- FreezeTech™ indicator changes color at 0°C
- Hi-Viz™ graduations with max freeze line
- Resealable sterile bags & eco-plastic racks

#### Eco-Plastic Racks



Autoclavable, recyclable,  
easily flattened after use



**BUY 2  
GET 1  
FREE**

## Solution Reservoirs

### Standard 12-Channel Width

**CLICK HERE  
to view these  
promos online**



50 mL - Shown with lid



25 mL



50 mL  
Divided



50 mL



100 mL

- Graduated inner walls
- Micro-Trough U-shaped bottoms
- Pouring spouts on all 4 corners
- RNase, DNase, & pyrogen free
- Pack of 5 protective lids with standard 25 and 50mL reservoirs

Cat No.	Description	Qty
21-415	15mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500 / case
21-416	15mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
21-417	50mL, Sterile, Hi-Viz™ Graduations, Eco-Plastic Racks of 25 in Resealable Bags	500 / case
21-418	50mL, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case

Description	25mL Standard*	50mL Standard	50mL Divided	100mL Standard	QTY
Non-Sterile Bulk Pack	28-134	28-109	28-160	28-163	100 / case
Sterile Individually Wrapped	28-134B	28-110	28-161	28-164	100 / case
Sterile Packs of 5	28-134C	28-111	28-162	28-165	200 / case
Lids - Sterile Packs of 5	28-135	28-111LID	N/A	N/A	200 / case

\*5 FREE Lids included with every purchase of Standard 25mL and Standard 50mL reservoirs.

## Poseidon™ Digital

### Electronic Pipette Controller



- Enhanced flow control – Max aspirate and dispense rates
- Independent aspiration and dispense speed settings
- Pressure sensitive triggers for even more control
- Improved stand - Locks in place with a click
- USB-C charge port – No searching for specific chargers



33-901D  
Digital Controller

Cat No.	Description	Qty
33-901D	Poseidon™ Digital Pipette Controller, Retractable Stand, USB-C Plug	1 each

**B2G1  
FREE**

**NEW!**

## Pipette Controllers

### Manual, Electronic, Long-Neck

#### Manual Controller 33-902

- No batteries required
- Autoclavable lower body
- Single-lever controls

#### Electronic Controller 33-901

- Retractable kickstand
- Contoured handle
- Available w/ US or EU plug

#### Long-Neck Controller 33-901XL

- Reduced arm strain
- Perfect for under the hood
- Available w/ US or EU plug



**B2G1  
FREE**

Cat No.	Description	Qty
<b>Other Poseidon Pipette Controllers</b>		
33-901	Poseidon™ Pipette Controller, Retractable Stand, US Type Plug	1 each
33-901XL	Poseidon™ Pipette Controller XL, Retractable Stand, US Type Plug	1 each
33-902	Poseidon™ Manual Pipette Controller, Single Lever Controller	1 each

To receive promotional packages, please reference **Promo Code RE2025** on your order. All promotions are valid while supplies last and expire **December 31, 2025**.  
All prices, specifications, and offers effective **April 1, 2025**. Subject to change without notice.

## BenchBin™ Biohazard Bin

Benchtop Wastebin with Spring Lid

- Safely dispose of tips, tubes, gloves, etc.
- Secures and covers hazardous waste
- Easy push open/close spring lid
- Compliant with OSHA standard 29 CFR 1910.1030
- Compact footprint – Minimize space on the bench



Cat. No.	Description	Qty
30-165	BenchBin™, 7×5×7", 1.5L Capacity, fits 30-162R Autoclave Bags	1 each
30-166	BenchBin 3L™, 6×5×8½", 3L Capacity, fits 30-166R Autoclave Bags	1 each
30-162R	Autoclave Bags, 8½×11" Printed Markings, fits 30-165 BenchBin™	100 / pack
30-166R	Autoclave Bags, 10×12½" Printed Markings, fits 30-166 BenchBin 3L™	100 / pack
30-165K	Starter Kit: 1 BenchBin™ 1.5L (30-165) and 400 Autoclave Bags (30-162R)	1 each
30-166K	Starter Kit: 1 BenchBin 3L™ (30-166) and 400 Autoclave Bags (30-166R)	1 each

## BowTie™ Biohazard Bags & Bins

Drawstring Bags and Hands-Free Bins



- Easy drawstring closure  
*No wire tie needed!*
- 12 or 30 gallon capacity bags
- Color-changing sterilization indicator
- Hands-free foot pedal lid on bins



Cat. No.	Description	Qty
<b>BowTie™ Biohazard Bags - Drawstring, Marking Area, Sterilization Indicator</b>		
30-210	27×35", Fits 30-211 and 30-219 Biohazard Bin	100 / pack
30-218	33×41", Fits 08-026 Biohazard Bin	100 / pack
<b>BowTie™ Biohazard Bins - Attached Lid and Hands-Free Foot Pedal</b>		
30-211	8 Gallon, 10×14×19", Fits 30-210 Biohazard Bag	1 Each
30-219	13 Gallon, 16×16×24", Fits 30-210 Biohazard Bag	1 Each
08-026	26 Gallon, 19×19×33", Fits 30-218 Biohazard Bag	1 Each

**NEW!**

## Olympus Plastics™ FACS Tubes

Tubes and Strainer Caps for Flow Cytometry

### Material

Tube: Polystyrene  
Cap Frame: HDPE  
Cap Mesh: Nylon

- Compatible with FACS™ (BD) instruments
- Strainer Caps: 20µm, 35µm, 70µm, 100µm mesh  
*Standard, non-strainer caps also available*
- 12×75mm (5mL) tubes in crystal clear polystyrene



Cat No.	Description	Qty
<b>Strainer and Standard Caps without Tubes</b>		
28-157	20µm, Red, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-155	35µm, Green, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-158	70µm, White, Strainer Cap, Opaque, Sterile, Bags of 25	500 / pack
28-159	100µm, Yellow, Strainer Cap, Translucent, Sterile, Bags of 25	500 / pack
28-152	Standard Cap, Non-Strainer, Non-Sterile, Bags of 125	1000 / case
<b>FlowTubes™ - FACS Tubes with Included Caps</b>		
28-154	5mL, 12×75mm, 35µm Strainer Cap, Sterile, Bags of 25	500 / case
28-156	5mL, 12×75mm, 35µm Strainer Cap, Sterile, Individual Wrap	500 / case
28-150	5mL, 12×75mm, Standard Cap, Sterile, Bags of 25	500 / case
<b>FlowTubes™ - FACS Tubes without Caps</b>		
28-151	5mL, 12×75mm, Without Cap, Sterile, Tray Packs of 250	1000 / case
28-153	5mL, 12×75mm, Without Cap, Non-Sterile, Bulk Econo-Pack	4000 / case





## GenClone® Cell Strainers

Pore Sizes Ranging from 5µm to 100µm

- 8 pore sizes for greater cell sorting precision
- Advanced design - faster flow through
- Included adapter for tubes smaller than 50mL
- Compatible with 6-well plates
- Interlocking and stackable
- 1.5mm feet for use in dishes & plates



5µm Strainer  
25-381 - Orange



10µm Strainer  
25-382 - Teal



20µm Strainer  
25-383 - Red



30µm Strainer  
25-384 - Purple



40µm Strainer  
25-375 - Blue



50µm Strainer  
25-385 - Grey



70µm Strainer  
25-376 - White



100µm Strainer  
25-377 - Yellow



Cat No.	Description	Qty
<b>GenClone® Advanced Cell Strainers</b>		
25-381	Cell Strainer, 5µm, Orange, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-382	Cell Strainer, 10µm, Teal, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-383	Cell Strainer, 20µm, Red, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-384	Cell Strainer, 30µm, Purple, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-375	Cell Strainer, 40µm, Blue, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-385	Cell Strainer, 50µm, Grey, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-376	Cell Strainer, 70µm, White, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
25-377	Cell Strainer, 100µm, Yellow, Sterile, Individually Wrapped, w/ Reducing Adapter	50/pk
<b>GenClone® Cell Strainer Accessories</b>		
25-278	Reducing Adapter, Sterile, Fits 5mL, 15mL, 12×75mm, and 17×100mm	25/pk
25-280	PolyPestle™ Double-Sided Pestle for Cell Strainers	50/pk
21-395	Centrifuge Tube, Attached Cap, Sterile, Self-Standing, Bags of 25	500/cs

NEW!

NEW!

NEW!

NEW!

NEW!

Reducing Adapter - 25-278



Allows coupling with 5mL, 15mL,  
12×75mm and 17×100mm Tubes

Individual Blister Packaging



Sterile; RNase, DNase, pyrogen  
and cytotoxin free




## Centrifuge Tubes

5mL to 500mL Capacities

- Perfect size for any application
- RNase, DNase, & pyrogen free
- Vivid graduations

Self-Standing

Self-Standing



Temperature

Min: -86°C

Max: 121°C

Material

Tube: Polypropylene

Cap: Polyethylene

Cat No.	Description	Qty
<b>5mL MacroTubes®</b>		
24-285SC	Separate Caps, Non-Sterile, Bulk Bagged	500 / case
24-285SCB	Attached Caps, Sterile, Bags of 100	500 / case
24-285SCS	Attached Caps, Sterile, Racks of 50	500 / case
24-285SBK	Attached Caps, Opaque Black, Sterile, Racks of 50	500 / case
24-285V	Attached Caps, Vivid Graduations, Sterile, Racks of 50	500 / case

### 10mL Centrifuge Tubes

21-412	Attached Cap, Sterile, Bags of 50	1000 / case
--------	-----------------------------------	-------------

### 15mL Centrifuge Tubes

21-408	Attached Cap, Sterile, Racks of 50	500 / case
21-408B	Attached Cap, Sterile, Bags of 50	500 / case
21-419	Attached Cap, Sterile, Eco-Plastic Racks of 25	500 / case
21-400	Attached Cap, Sterile, Individually wrapped	300 / case

Cat No.	Description	Qty
<b>25mL Centrifuge Tubes</b>		
24-295	Attached Cap, Sterile, Racks of 25	200 / pack
24-296	Attached Cap, Sterile, Bags of 25	200 / pack

### 30mL Centrifuge Tubes

21-402	Screw Cap, Self-Standing, Sterile, Bags of 25	500 / case
--------	---	------------

### 50mL Centrifuge Tubes

21-409	Attached Cap, Sterile, Racks of 25	500 / case
21-409B	Attached Cap, Sterile, Bags of 25	500 / case
21-421	Separate Screw Cap, Non-Sterile, Bulk Bag	500 / case
21-420	Attached Cap, Sterile, Eco-Plastic Racks of 25	500 / case
21-401	Attached Cap, Sterile, Individually Wrapped	300 / case
21-395	Attached Cap, Sterile, Self-Standing, Bags of 25	500 / case


Cat No.	Description	Qty
<b>Large Volume Centrifuge Tubes</b>		
21-422	175mL, Sterile PP, Conical Bottom, Bags of 6	120 / case
21-423	225mL, Sterile, PP, Conical Bottom, Bags of 5	90 / case
21-404	250mL, Sterile, PP, Conical Bottom, Bags of 5	70 / case
08-068	500mL, Sterile, PP, Conical Bottom, Bags of 5	40 / case

Cat No.	Description	Qty
<b>Centrifuge Tube Racks</b>		
27-470	4×50/25mL, 12×15/10/5mL, 32×1.5/0.5mL Tubes	5 / pack
27-117	Decanting, 50×5mL, Silicone Grip Wells	1 each
27-386	30 Slots for 10/15mL, 20 Slots for 25/50mL	5 / pack
27-471	6 Slots for 175, 225, or 250mL Tubes	5 / pack
24-323	4 Slots for 500mL Tubes	1 each

## Disposable Cuvettes

### Spectrophotometry Cuvettes

- Fits most spectrophotometers and photometers
- Polystyrene for use in visible range (340-800nm)
- Protective packaging – dust-proof and non-scratch
- Absorption reproducibility error of less than  $\pm 1\%$



Cat No.	Description	Qty	
21-190	Cuvette, PS, 1.5mL semi-micro, round opening, 10mm light path	1000 / pack	NEW!
21-191	Cuvette, PS, 2.5mL semi-micro, round opening, 10mm light path	1000 / pack	NEW!
21-192	Cuvette, PS, 1.5mL semi-micro, square opening, 10mm light path	1000 / pack	NEW!
21-193	Cuvette, PS, 4.5mL macro, square opening, 10mm light path	1000 / pack	NEW!
21-194	Cuvette Plug Cap, Round Closure, fits S4100-1 and S4100-2 cuvettes	500 / pack	NEW!
21-195	Cuvette Plug Cap, Square Closure, fits S4100-3 and S4100-4	500 / pack	NEW!

## Olympus Reservoirs

### Reduced 8-Channel Width

- Specifically designed for 8-channel pipettes
- 75mm inner chamber
- Conserve valuable reagents




Description	5mL 8-Channel	10mL 8-Channel	25mL 8-Channel	QTY
Non-Sterile Bulk Pack	28-112	28-115	28-131*	300 / case*
Sterile Individually Wrapped	28-113	28-116	28-132	100 / case
Sterile Packs of 5	28-114	28-117	28-133	200 / case

\*28-131 has quantity of 100/pack

## ClearSeal™ Screw Cap Microtubes

- Leak-proof O-Ring seal
- Ultra-clear polypropylene
- Resealable sterile bags of 50
- RNase, DNase, and Pyrogen free
- White printed graduations
- Conical or standing bottoms




Temperature	Material	Max RCF
Min: -86°C Max: 121°C	Polypropylene	22,000xg

Cat No.	Description	Qty
24-292	0.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000 / case
24-292B	0.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap,	1000 / case
24-297	1.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000 / case
24-297B	1.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000 / case
24-298	2.0mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap, Graduated	1000 / case
24-298B	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000 / case
24-299	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Loop Screw Cap, Graduated	1000 / case

## TotalPure™ Tubes

### Metal-Free Centrifuge Tubes

- Ultra-Pure polypropylene  
- Certified free of trace metals
- Hi-Viz™ graduations  
- w/ Max Freeze Line
- 20,000xg Max RCF
- Sterile - resealable bags or individually wrapped




Sterile Packaging:  
Individual or Bags of 25

Cat No.	Description	Qty
NEW! 21-425	15mL, Metal Free, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
NEW! 21-426	50mL, Metal Free, Sterile, Hi-Viz™ Graduations, Resealable Bags of 25	500 / case
NEW! 21-427	15mL, Metal Free, Sterile, Hi-Viz™ Graduations, Individually Wrapped	500 / case
NEW! 21-428	50mL, Metal Free, Sterile, Hi-Viz™ Graduations, Individually Wrapped	500 / case

## OptiWell™ Consumables

### For KingFisher™ Platforms



### Deep Well Plates

Printed Alphanumerics  
Virgin, medical grade PP  
24 or 96 square well models

### Sealing Mats

Leak-proof silicone sealing  
Printed Alphanumerics  
Minimize sample loss

### Tip Combs

Medical-grade PP  
KingFisher™ Flex compatible  
8, 24 or 96 tip models

Cat No.	Description	Qty
<b>Deep Well Plate - Square Well, V-Bottom, PP, Printed Alpha-numerics</b>		
30-401	Deep Well Plate, 24 Square Well, 10mL, 41.9mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
30-415	Elution Plate, 96 Square Well, 0.5mL, 15.1mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
30-425	Deep Well Plate, 96 Square Well, 1.0mL, 24.4mm Height, V-Bottom, Sterile, for KingFisher™	100 / case
30-445	Deep Well Plate, 96 Square Well, 2.2mL, 44.2mm Height, V-Bottom, Sterile, for KingFisher™	50 / case
<b>Sealing Mat - Square Well, Silicone, Printed Alpha-numeric</b>		
30-467	Sealing Mat, 24 Square Well, Silicone, Sterile	50 / case
30-455	Sealing Mat, 96 Square Well, Ø8.4mm, "+" Opening, Alphanumerics, Silicone, Sterile	50 / case
30-457	Sealing Mat, 96 Square Well, Ø8.4mm, "*" Opening, Alphanumerics, Silicone, Sterile	50 / case
<b>Tip Comb - V-Bottom, PP</b>		
30-473	Tip Comb, 8 Tip, V-Bottom, Sterile, PP, for KingFisher™	480 / case
30-475	Tip Comb, 24 Tip, V-Bottom, Sterile, PP, for KingFisher™	50 / case
30-477	Tip Comb, 96 Tip, V-Bottom, Sterile, PP, for KingFisher™	100 / case

Visit [geneseesci.com](http://geneseesci.com) for non-sterile models and dozens of other plate/mat configurations.

KingFisher™ is a registered trademark of Thermo Fisher Scientific

## PCR Sealing Films & Foils

### Sterile Seals for PCR, qPCR, ELISA, EIA

- All types available sterile or non-sterile
- Evaporation rates < 5%
- 10 inner bags of 10, 100 per pack
- RNase, DNase, DNA, & PCR inhibitor free



Cat. No.	Description	Qty
NEW! 12-205	12-205S PCR Sealing Foil, Aluminum Foil, Adhesive, Pierceable	100 / pack
NEW! 12-207	12-207S Sealing Film, Transparent, Easy Peel Low Adhesive, for ELISA & EIA	100 / pack
NEW! 12-209	12-209S PCR Sealing Film, Transparent, Adhesive	100 / pack
NEW! 12-210	12-210S qPCR Sealing Film, Optically Clear, Standard Adhesive	100 / pack
NEW! 12-212	12-212S qPCR Sealing Film, Optically Clear, Non-Tacky Pressure Activated Adhesive	100 / pack