



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Haier - DW-86L959BPT(60Hz) : DW-86L959BPT(60Hz)

Specifications

ENERGY STAR Unique ID:	2384620
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-86L959BPT(60Hz)
Model Number:	DW-86L959BPT(60Hz)
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	33.9
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.33
Peak Temperature Variance (°C):	16.8
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Right,Swinging
Number of Drawers or Shelves:	3
Number of Internal Ports and/or Access Holes:	4
Depth (in):	39.29
Height:	77.95
Width:	51.02
Defrost Type:	Manual
Refrigerant with GWP:	R-170 (ethane) (GWP:6 Lower GWP), R-290 (GWP:3 Lower GWP)
Installed Accessories:	N/A
Date Available on Market:	2021-09-17
Date Certified:	2021-09-23
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Additional Model Information

„100-230V,60Hz

UPC Codes

Captured On:
01/19/2022