



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Haier - DW-86L728J : DW-86L728J

Specifications

ENERGY STAR Unique ID:	2371855
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-86L728J
Model Number:	DW-86L728J
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	25.42
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.38
Peak Temperature Variance (°C):	14.2
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging,Left
Number of Drawers or Shelves:	3
Number of Internal Ports and/or Access Holes:	4
Depth (in):	38.58
Height:	77.95
Width:	41.3
Defrost Type:	Manual
Refrigerant with GWP:	R-170 (ethane) (GWP:6 Lower GWP), R-290 (GWP:3 Lower GWP)
Installed Accessories:	N/A
Date Available on Market:	2021-01-10
Date Certified:	2021-01-12
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Additional Model Information

UPC Codes

Captured On:
01/19/2022