



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Haier - DW-86L728J(220V) : DW-86L728J(220V)

Specifications

ENERGY STAR Unique ID:	2321351
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-86L728J(220V)
Model Number:	DW-86L728J(220V)
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	25.32
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.42
Peak Temperature Variance (°C):	21.18
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Right,Swinging
Number of Drawers or Shelves:	4
Number of Internal Ports and/or Access Holes:	2
Depth (in):	37.25
Height:	78.0
Width:	37.5
Defrost Type:	Manual
Refrigerant with GWP:	R-170 (ethane) (GWP:6 Lower GWP), R-290 (GWP:3 Lower GWP)
Installed Accessories:	Chart recorder
Date Available on Market:	2014-06-30
Date Certified:	2018-06-15
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Additional Model Information

„220-240V, 50H
Z

UPC Codes

Captured On:

01/19/2022