



ENERGY STAR CERTIFIED

# Lab Grade Refrigerators and Freezers

## Haier - DW-86L828JA : DW-86L828JA

### Specifications

ENERGY STAR Unique ID:	2376352
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-86L828JA
Model Number:	DW-86L828JA
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	29.09
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.36
Peak Temperature Variance (°C):	15.6
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging,Left
Number of Drawers or Shelves:	3
Number of Internal Ports and/or Access Holes:	4
Depth (in):	38.58
Height:	77.95
Width:	45.08
Defrost Type:	Manual
Refrigerant with GWP:	R-170 (ethane) (GWP:6   Lower GWP), R-290 (GWP:3   Lower GWP)
Installed Accessories:	N/A
Date Available on Market:	2021-04-20
Date Certified:	2021-04-26
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

### Additional Model Information

,,208-230V~,60  
Hz

UPC Codes

Captured On:

01/19/2022