



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Haier - DW-30L818BP (60Hz) : DW-30L818BP (60Hz)

Specifications

ENERGY STAR Unique ID:	2361986
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-30L818BP (60Hz)
Model Number:	DW-30L818BP (60Hz)
Product Type:	High Performance Freezer
Total Volume (cubic feet):	28.9
Lab Grade Refrigerator or Freezer Energy Consumption (kWh/day):	4.2
Peak Temperature Variance (°C):	9.6
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging,Left
Number of Drawers or Shelves:	5
Number of Internal Ports and/or Access Holes:	3
Depth (in):	29.72
Height:	77.95
Width:	29.53
Defrost Type:	Manual
Refrigerant with GWP:	R-290 (GWP:3 Lower GWP)
Installed Accessories:	N/A
Date Available on Market:	2020-06-22
Date Certified:	2020-06-23
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Additional Model Information

„100-230V 60Hz

UPC Codes

Captured On:
01/19/2022