

HandyStep®

Repeating Pipettes



BRAND® HandyStep® electronic

The HandyStep® electronic motorized repeating pipette combines broad tip compatibility and user comfort with the convenience and efficiency of an electronic pipette.

- **Works with All Standard Tips:** The HandyStep® electronic is compatible with BRAND PD-Tip™ syringe tips, Eppendorf® Combitips advanced®, and more.
- **Recognizes Tips Automatically:** The pipette automatically recognizes the size of encoded BRAND PD-Tip™ syringe tips upon tip insertion and displays the tip's volume, default dispensing volume, and default number of dispensing steps.
- **Reduces Repetitive Motion by up to 97%:** The instrument features a unique Auto-Dispense operating mode in which the pipette “learns” the operator’s dispensing rhythm and dispenses at that rhythm as long as the control button is depressed – no programming needed!

Use the tips you want

The only electronic repetitive pipette compatible with all standard syringe tips, the HandyStep® electronic accepts:

- BRAND PD-Tip™ syringe tips
- Rainin® Encode™ syringe tips
- Eppendorf® Combitips® advanced
- Fisherbrand® Dispenser Tips
- VWR® Combi-Syringes
- Gilson® Repet-Tips

Description	Cat. No.	2017 List Price
HandyStep® electronic with charging stand & NiMH battery, 110V, 50-60Hz	705002	\$1,019.00
HandyStep® S Repeating Pipette	705110	378.00

BRAND® HandyStep® S

The HandyStep® S repeating pipette provides unprecedented versatility in a purely mechanical repeating pipette. In conjunction with BRAND PD-Tip™ syringe tips, the HandyStep® S offers nearly sixty different volume settings with accuracy, precision and dependability.

- **Comfortable:** The slim design of the HandyStep® S is lightweight, well-balanced, and provides a convenient control layout for smooth dispensing and one handed volume selection.
- **Reliable:** Robust construction, inspired design and innovative polymers deliver a dependable instrument.
- **Ready-To-Go:** A purely mechanical design with no batteries to charge or replace means the HandyStep® S is ready when you are. A convenient integrated volume table provides a quick reference for dispensed volumes.
- **Flexible:** Ten sizes of PD-Tip™ syringe tips allow the HandyStep® S to dispense 59 different volumes. Open design allows use of third party tips—no more being locked into a sole source!
- **Contact-Free:** A tip ejection system allows the user to remove tips without handling.

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726
www.brandtech.com

PD-Tip™ Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette.

BRAND encoded PD-Tip™ syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** PD-Tip™ syringe tips are compatible with repeating pipettes using standard tips. PD-Tip™ syringe tips are encoded for automatic recognition by the BRAND HandyStep® electronic repeating pipette, the Rainin® AutoRep™ E pipette, and the Gilson® Repetman®.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip™ syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **BIO-CERT® PD-Tip™ syringe tips** are sterile, individually wrapped and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of BRAND PD-Tip™ syringe tips.



PD-Tip™ Syringe Tips, Syringe Tips, BIO-CERT® sterile, individually wrapped

Capacity, mL	Pack of	Cat. No.	2017 List Price
0.1mL	100	702683	\$167.00
0.5mL	100	702684	160.00
1mL	100	702685	160.00
1.25mL	100	702686	160.00
2.5mL	100	702688	160.00
5mL	100	702690	160.00
10mL	100	702691	160.00
12.5mL	100	702692	160.00
25mL, (one adapter included)	25	702694	47.00
50mL, (one adapter included)	25	702696	48.00
Adapter for 25mL and 50mL PD-Tip™ BIO-CERT® sterile, PP, autoclavable	5	702699	74.00

*The Eppendorf® Repeater® Plus, Repeater® M4, Repeater® Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf® Combitips® advanced.

PD-Tip™ Syringe Tips, non-sterile

Capacity, mL	Pack of	Cat. No.	2017 List Price
0.1mL	100	702402	\$125.00
0.5mL	100	702370	117.00
1mL	100	702406	117.00
1.25mL	100	702372	117.00
2.5mL	100	702374	117.00
5mL	100	702376	117.00
10mL	100	702407	117.00
12.5mL	100	702378	117.00
25mL, (one adapter included)	50	702380	73.00
50mL, (one adapter included)	25	702382	42.00
PD-Tip™ Syringe Tips, non-sterile, assorted pack 0.5, 1, 1.25, 2.5, 5, 10 and 12.5mL sizes	20/size	702368	158.00
Adapter for 25mL and 50mL PD-Tip™ Syringe Tips non-sterile, PP, autoclavable	10	702398	94.00

HandyStep_BT_Fly-2017



Lab Rats Trust BrandTech!

Toll Free: 888-522-2726 www.brandtech.com

BRANDTECH®
SCIENTIFIC, INC.