



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Control Normal Sheep IgG

CATALOG NUMBER: 5091

Specifications

Host Species	Sheep
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC, WB

Properties

Purification	Normal sheep IgG is obtained from serum of sheep that have not been immunized. This preparation can be used as a control reagent for immunoassays using sheep polyclonal antibodies. The IgG was purified by DEAE-column chromatography.
Isotype	IgG
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Stable at 4°C for three months, store at -20°C for up to one year.

Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).