

TCTN3 Antibody


CATALOG NUMBER: 5837

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	TCTN3 antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of human TCTN3. The immunogen is located within amino acids 540 - 590 of TCTN3.
Conjugate	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	TCTN3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	TCTN3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	TCTN3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.




TCTN3 Antibody 1

Western blot analysis of TCTN3 in HeLa cell lysate with TCTN3 antibody at (A) 1 and (B) 2µg/mL.



TCTN3 Antibody 2

Immunocytochemistry of TCTN3 in HeLa cells with TCTN3 antibody at 4 µg/mL.



TCTN3 Antibody 3

Immunofluorescence of TCTN3 in HeLa cells with TCTN3 antibody at 20 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).