



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## Control Normal Monkey IgG

CATALOG NUMBER: 5891

### Specifications

Host Species	Monkey
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC, WB

### Properties

Purification	Normal monkey IgG is obtained from serum of monkeys that have not been immunized. This preparation can be used as a control reagent for immunoassays. The IgG was purified by Protein A-column chromatography.
Isotype	IgG
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Stable at 4°C for three months, store at -20°C for up to one year.

### Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).