



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Pancreas Lysate

CATALOG NUMBER: 1414

Specifications

Species	Mouse
Diagnosis	Normal
Tested Applications	WB

Properties

Physical State	Liquid
Buffer	Pancreas lysate is supplied in SDS sample buffer containing 50 mM DTT.
Concentration	2 mg/mL
Storage Conditions	Store at 2-8°C for continuous use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt.

Disclaimer

Disclaimer	This product is intended for laboratory research purposes only and should be used by qualified personnel only. ProSci lysates are not intended for use in humans. ProSci is not liable for damages or injuries resulting from receipt and/or use of ProSci lysates.
------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).