

TCTEX1D3 Antibody


CATALOG NUMBER: 6955

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	TCTEX1D3 antibody was raised against a 19 amino acid synthetic peptide near the amino terminus of human TCTEX1D3. The immunogen is located within the first 50 amino acids of TCTEX1D3.
Conjugate	Unconjugated
Tested Applications	ELISA, IF, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	TCTEX1D3 is predicted to not cross-react with other members of the TCTEX-type family
Predicted Molecular Weight	22 kDa


Properties

Purification	TCTEX1D3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	TCTEX1D3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	TCTEX1D3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



TCTEX1D3 Antibody 1

Western blot analysis of TCTEX1D3 in EL4 cell lysate with TCTEX1D3 antibody at 1 µg/mL in (A) the absence and (B) the presence of blocking peptide.



TCTEX1D3 Antibody 2

Immunofluorescence of TCTEX1D3 in human testis tissue with TCTEX1D3 antibody at 20 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).