

## TDGF1 Antibody


CATALOG NUMBER: 7099

### Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	Rabbit polyclonal TDGF1 antibody was raised against a 17 amino acid peptide near the amino terminus of human TDGF1. The immunogen is located within amino acids 40 - 90 of TDGF1.
Conjugate	Unconjugated
Tested Applications	ELISA, IF, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	21 kDa


### Properties

Purification	TDGF1 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	TDGF1 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	TDGF1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



TDGF1 Antibody 1

Western blot analysis of TDGF1 in 293 cell lysate with TDGF1 antibody at 1 µg/mL.



TDGF1 Antibody 2

Immunofluorescence of TDGF1 in 293 cells with TDGF1 antibody at 20 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).