

VTCN1 Antibody


CATALOG NUMBER: 7189

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse
Homology	Predicted species reactivity based on immunogen sequence: Rat: (83%)
Immunogen	Rabbit polyclonal VTCN1 antibody was raised against an 18 amino acid peptide near the center of human VTCN1. The immunogen is located within amino acids 70 - 120 of VTCN1.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	Predicted: 31 kDa

Properties

Purification	VTCN1 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	VTCN1 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	VTCN1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



VTCN1 Antibody 1

Western blot analysis of VTCN1 in EL4 cell lysate with VTCN1 antibody at 1 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).