

BTLA Antibody


CATALOG NUMBER: 7473

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	BTLA antibody was raised against a 17 amino acid peptide near the center of human BTLA. The immunogen is located within amino acids 170 - 220 of BTLA .
Conjugate	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	BTLA antibody is human specific. At least two isoforms of BTLA are known to exist.
Predicted Molecular Weight	Predicted: 32 kDa Observed: 40kDa


Properties

Purification	BTLA Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	BTLA Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	BTLA antibody can be stored at 4°C for three months and -20°C, stable for up to one year.




BTLA Antibody 1

Western blot analysis of BTLA in Jurkat cell lysate at (A) 1 and (B) 2 $\mu\text{g}/\text{mL}$.



BTLA Antibody 2

Immunohistochemistry of BTLA in human small intestine tissue with BTLA antibody at 2.5 $\mu\text{g}/\text{mL}$.



BTLA Antibody 3

Immunofluorescence of BTLA in human small intestine tissue with BTLA antibody at 20 $\mu\text{g}/\text{mL}$.

Disclaimer

Disclaimer

Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).