



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Andes Virus Nucleocapsid Antibody

CATALOG NUMBER: 7679

Specifications

Host Species	Rabbit
Species Reactivity	Virus
Immunogen	Andes virus nucleocapsid antibody was raised against a 19 amino acid peptide near the amino terminus of the Andes virus nucleocapsid. The immunogen is located within amino acids 130 - 180 of Andes Virus Nucleocapsid.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	N/A

Properties

Purification	Andes virus nucleocapsid antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	Andes virus nucleocapsid antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Andes virus nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).