



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Mouse Anti-Pig IgG Antibody [1G5H7] (biotin)

CATALOG NUMBER: 1608

Specifications

Host Species	Mouse
Species Reactivity	Pig
Immunogen	Pig IgG
Conjugate	Biotin
Tested Applications	ELISA, IHC, WB
User Note	Prepare working dilution fresh each time.
Specificity	No cross reactivity with immunoglobulins from other species.

Properties

Purification	The purified antibody-Biotin conjugate is separated from non-conjugated biotin by dialysis.
Isotype	IgG1
Physical State	Liquid
Buffer	PBS pH 7.2, 1% BSA, 0.02% Sodium Azide
Concentration	1 mg/mL
Storage Conditions	Store at 2-8°C. The product is stable for about 4-6 weeks at 4°C as an undiluted liquid. For extended storage, aliquot and freeze the product at -70°C or below. Avoid repeated freezing and thawing. Or alternatively, add an equal volume of glycerol for a final concentration of 50% and store at -20°C.

Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).