



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Control Normal Llama IgG1

CATALOG NUMBER: 8399

Specifications

Host Species	Llama
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC, WB

Properties

Purification	Normal Llama IgG1 is obtained from serum of Llamas that have not been immunized. This preparation can be used as a control reagent for immunoassays using Llama polyclonal antibodies. The IgG1 was purified by Protein G-column chromatography.
Isotype	IgG1
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Stable at 4°C for three months, store at -20°C for up to one year.

Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).