



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

Control Normal Llama IgG (Total)

CATALOG NUMBER: 8411

Specifications

| | |
|---------------------|-------------------|
| Host Species | Llama |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | ELISA, IHC, WB |

Properties

| | |
|--------------------|--|
| Purification | Normal Llama IgG is obtained from serum of Llamas that have not been immunized. This preparation can be used as a control reagent for immunoassays using Llama polyclonal antibodies. The IgG was purified by Protein A-column chromatography. |
| Isotype | IgG |
| Physical State | Liquid |
| Buffer | PBS containing 0.02% sodium azide. |
| Concentration | 1 mg/mL |
| Storage Conditions | Stable at 4°C for three months, store at -20°C for up to one year. |

Disclaimer

| | |
|------------|--|
| Disclaimer | This product is for research use only. |
|------------|--|

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).