

UEV1A Antibody


CATALOG NUMBER: 3375

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Homology	Predicted species reactivity based on immunogen sequence: Bovine: (100%), Chicken: (100%)
Immunogen	UEV1A antibody was raised against a peptide corresponding to 15 amino acids near the C-terminus of human UEV1A. The immunogen is located within the last 50 amino acids of UEV1A.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	Anti-UEV1A may also recognize other isoforms of UEV1A.

Properties

Purification	UEV1A Antibody is Ion exchange chromatography purified.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	UEV1A Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	UEV1A antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



UEV1A Antibody 1

Western blot analysis of UEV1A in Jurkat cell lysates with UEV1A antibody at (A) 1, (B) 2, and (C) 4 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).