



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

FOXO3 Antibody


CATALOG NUMBER: 5813

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	FOXO3 antibody was raised against an 18 amino acid synthetic peptide near the amino terminus of human FOXO3. The immunogen is located within amino acids 40 - 90 of FOXO3.
Conjugate	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	FOXO3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	FOXO3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	FOXO3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



FOXO3 Antibody 1

Western blot analysis of FOXO3 in A-20 cell lysate with FOXO3 antibody at (A) 0.5 and (B) 1 µg/mL.



FOXO3 Antibody 2

Immunocytochemistry of FOXO3 in A-20 cells with FOXO3 antibody at 4 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).