

ARHGAP39 Antibody


CATALOG NUMBER: 6445

Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	ARHGAP39 antibody was raised against an 18 amino acid synthetic peptide near the carboxy terminus of human ARHGAP39. The immunogen is located within amino acids 880 - 930 of ARHGAP39.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	ARHGAP39 antibody is predicted to not cross-react with other RhoGAP protein family members.

Properties

Purification	ARHGAP39 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	ARHGAP39 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	ARHGAP39 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



ARHGAP39 Antibody 1

Western blot analysis of ARHGAP39 in A20 cell lysate with ARHGAP39 antibody at (A) 1 and (B) 2 µg/mL

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).