

BCL3 Antibody


CATALOG NUMBER: 7183

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	Rabbit polyclonal BCL3 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human BCL3. The immunogen is located within amino acids 360 - 410 of BCL3.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	Predicted: 50 kDa

Properties

Purification	BCL3 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	BCL3 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	BCL3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



BCL3 Antibody 1

Western blot analysis of BCL3 in 293 cell lysate with BCL3 antibody at 1 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).