

## Heartland Virus Glycoprotein 2 Antibody

CATALOG NUMBER: 7435

### Specifications

Host Species	Rabbit
Species Reactivity	Virus
Immunogen	Heartland virus glycoprotein 2 antibody was raised against a 19 amino acid peptide near the center of the Heartland virus glycoprotein 2. The immunogen is located within amino acids 860 - 910 of Heartland Virus Glycoprotein 2.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	Predicted: 55 kDa

### Properties

Purification	Heartland Virus Glycoprotein 2 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	Heartland virus glycoprotein 2 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Heartland virus glycoprotein 2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

### Disclaimer

#### Disclaimer

Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only. For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).