

Vinculin Antibody (biotin)


CATALOG NUMBER: 7807-biotin

Specifications

Host Species	Rabbit
Species Reactivity	Chicken, Human, Mouse, Rabbit, Rat
Immunogen	Biotin-Vinculin antibody was raised against a 23 amino acid peptide near the amino terminus of human Vinculin.
Conjugate	Biotin
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	Biotin-Vinculin antibody is human, mouse, rat, rabbit and chicken reactive. At least three isoforms of Vinculin are known to exist.
Predicted Molecular Weight	Predicted: 117 kDa Observed: 117 kDa

Properties

Purification	Biotin-Vinculin antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	Biotin-Vinculin Antibody is supplied in PBS containing 1% BSA and 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Biotin-Vinculin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



Vinculin Antibody (biotin) 1

Western blot analysis of Vinculin in 293, A431, A549, Daudi, HeLa, HepG2, K562, 3T3, Raji, U937 cell lysate and mouse brain, rat brain, rabbit brain, mouse lung, rat lung, mouse liver, rat liver, r...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).