

## Avian Influenza A H7N9 Neuraminidase Antibody


CATALOG NUMBER: 7831

### Specifications

Host Species	Rabbit
Species Reactivity	Virus
Immunogen	Avian Influenza A H7N9 Neuraminidase antibody was raised against a synthetic peptide corresponding to 11 amino acids near the amino terminus of the H7N9 [Influenza A virus (A/blue-winged teal/Ohio/566/2006(H7N9))] Neuraminidase protein. The immunogen is located within amino acids 50 - 100 of Avian Influenza A H7N9 Neuraminidase.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	virus

### Properties

Purification	Avian Influenza A H7N9 Neuraminidase antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	Avian Influenza A H7N9 Neuraminidase Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Avian Influenza A H7N9 Neuraminidase antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



Avian Influenza A H7N9  
Neuraminidase Antibody 1

H7N9 Neuraminidase antibody (Cat. No. 7831 at 1 µg/ml) detects 10 ng of H7N9 [Influenza A virus (A/blue-winged teal/Ohio/566/2006(H7N9))] Neuraminidase peptide, and not 10 ng of H7N9 [Influenza A v...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).