

SARS-CoV-2 (COVID-19) NSP2 Antibody


CATALOG NUMBER: 9171

Specifications

Host Species	Rabbit
Species Reactivity	Virus
Homology	Predicted reactivity based on immunogen sequence: SARS-CoV NSP2 protein: (79%)
Immunogen	Anti-SARS-CoV-2 (COVID-19) NSP2 antibody (9171) was raised against a peptide corresponding to 14 amino acids near the center of SARS-CoV-2 (COVID-19) NSP2 protein. The immunogen is located within 90-130 amino acids of SARS-CoV-2 (COVID-19) NSP2 protein.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	SARS-CoV-2 (COVID-19) NSP2 Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	SARS-CoV-2 (COVID-19) NSP2 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	SARS-CoV-2 (COVID-19) NSP2 antibody can be stored at 4 ° C for three months and -20 ° C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



**SARS-CoV-2 (COVID-19) NSP2
Antibody 1**

Figure 1 ELISA Validation with SARS-CoV-2 (COVID-19) NSP2 Protein
Antibodies: SARS-CoV-2 (COVID-19) NSP2 Antibody, 9171 A direct ELISA was performed using SARS-CoV-2 NSP2 rec...

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).